VETLIS® Brucella iELISA Caninos

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VETLIS®
Brucella iELISA Canine (serum)

Presentations / Components

Components	Description	15-VL03-2P 2 microplates	15-VL03-5P 5 microplates
Microplate	96-well microtitre plate (12 x 8 strip wells with strip holder) coated with the antigen, sealed and stored in dry	2 microplates (192 tests)	5 microplates (480 tests)
Conjugate 10X IgG	10X concentrated (horseradish peroxidase conjugated anti-canine IgG antibodies)	1 x 2.5 ml	1 x 6.5 ml
Substrate-chromogen solution	Ready to use [3,3',5,5'-Tetramethyl-bencidine (H ₂ O ₂)] – STORE IN THE DARK	1 x 25 ml	1 x 60 ml
Stop solution	Ready to use (contains HCl 1%) - CORROSIVE	1 x 25 ml	1 x 60 ml
Wash solution	10X concentrated	1 x 120 ml, 10X (c. s. p. 1.2 liters)	1 x 120 y 1 x 60 ml, 10X (c. s. p. 1.8 liters)
Sample diluent	Ready to use - CONTAINS PRE-SERVATIVE [Thimerosal (C ₉ H ₉ HgNaO ₂ S) 0.01%]	1 x 120 ml	2 x 120 y 1 x 60 ml
Positive control	Serum - CONTAINS PRESERVATIVE [Thimerosal (C ₉ H ₉ HgNaO ₂ S) 0.01%]	1 x 0.1 ml	2 x 0.1 ml
Negative control	Serum - CONTAINS PRESERVATIVE [Thimerosal (C ₉ H ₉ HgNaO ₂ S) 0.01%]	1 x 0.1 ml	2 x 0.1 ml
Polyethylene bag	With ziplock and reusable	1	1
Instruction manual		1	1

Name and application

VETLIS®

Brucella iELISA Canines (serum)

Indirect ELISA kit for the detection of anti *Brucella canis* antibodies in canine serum.

The kit VETLIS® *Brucella* iELISA Canine is an ELISA (Enzyme Linked Immunosorbent Assay) kit for the indirect detection of specific antibodies against the rough lipopolysaccharide (rLPS) of *Brucella canis* in canine serum.

Introduction

Brucellosis is a highly contagious zoonotic disease caused by Gram-negative bacteria of the genus *Brucella*. The species *Brucella canis* is the causative agent of brucellosis in dogs and also a human pathogen. *B. canis* causes abortions and stillbirths in pregnant bitches and epididymitis, orchitis and prostatitis in males. These clinical manifestations can lead to sterility, a factor that produces important economic losses in dog breeding kennels. Dogs that have been neutered do not have reproductive signs, but they occasionally develop other conditions such as ocular disease and discospondylitis. Among dogs the infection is transmitted through contaminated aborted fetuses, urine, milk, vaginal secretions and semen. Additionally, *B. canis* can spread through fomites. It has been demonstrated that uninfected dogs living with infected animals of the same sex can acquire the infection within 6 months.

The control of canine brucellosis depends primarily on the detection and separation of infected animals. The confirmatory diagnosis of brucellosis is performed directly by isolating the organism from blood and/or semen/vaginal swab culture. However, due to the slow growth of *Brucella* in primary cultures (up to 7 days), the risk involved in manipulating the bacteria and low sensitivity of its iso-

lation, determined that the isolation of *Brucella* for diagnosis is not always feasible and effective. For this reason, the laboratory diagnosis of canine brucellosis is mainly based on the detection of specific antibodies against the infectious agent in serum samples

Principles of the technique

The VETLIS® *Brucella* iELISA Canine is an enzyme immunoassay based on the detection of IgG antibodies directed against the *Brucella canis* rough lipopolysaccharide (rLPS).

In this assay, serum samples are exposed to the wells coated with the *Brucella canis* antigen (rLPS). If the sample contains anti rLPS antibodies, they bind to the antigen in the well. Upon addition of the horseradish peroxidase (HRP) conjugated antibody, a complex is formed with the canine IgG antibodies attached to the antigen in the well. Free antibodies are removed by washing before adding the substrate-chromogen solution. The blue color that appears is the product of the reaction of the substrate-chromogen with the conjugate in those samples that are positive. The color intensity is directly proportional to the amount of antibodies present in the sample reacting with the antigen. The higher the antibody titer present in the sample, the more intense color develops in the test wells. The reaction is stopped by adding the stop solution, upon which a color change to yellow is observed. The reading of the results is performed in a microplate spectrophotometer measuring the absorbance (Abs) at 450 nm.

Kit components

- Microplate
96-well microtitre plate (12 x 8 strip wells with strip holder) coated with the antigen, sealed and stored dry.
- Wash solution
10X concentrated.
- Conjugate 10X IgG

10X concentrated (horseradish peroxidase conjugated anti-canine IgG antibodies).

-Substrate-chromogen solution

Ready to use [3,3',5,5'-Tetramethylbenzidine in substrate buffer containing hydrogen peroxide (H_2O_2)]- STORE IN THE DARK.

-Stop solution

Ready to use (contains HCl 1%) - CORROSIVE.

-Sample diluent

Ready to use – CONTAINS PRESERVATIVE [Thimerosal ($C_9H_9HgNaO_2S$) 0.01%]. Shake well before use.

-Canine Positive Control Serum

Serum – CONTAINS PRESERVATIVE [Thimerosal ($C_9H_9HgNaO_2S$) 0.01%].

-Canine Negative Control Serum

Serum – CONTAINS PRESERVATIVE [Thimerosal ($C_9H_9HgNaO_2S$) 0.01%].

-Polyethylene bag

With ziplock and reusable.

-Instruction manual

Storage and expiration

Store: between 2 and 8°C.

Transport temperature: 4 to 15°C for 72 hours. Proper storage of the kit and all reagents will allow the use until the expiration date.

Expiration: 12 months. Do not use reagents from different lots or after the expiration date.

Materials needed but not provided

-High precision multichannel pipette and micropipettes (20, 100 and 1000 µl).

-Containers for multichannel pipette.

-Disposable pipette tips.

-Microplate reader for absorbance determination at 450 nm.

-Orbital rotator.

-Timer.

-Tubes for the dilutions of the controls and samples.

-1 liter container for the preparation of the wash solution.

-Distilled (Type II) or ultrapure (Type I) or any similar high quality water.

-Microplate covers (lid, aluminium foil or adhesive).

Precautions

1. Carefully read and strictly follow all instructions detailed in the Instruction manual included in the kit.

2. Store the kit and all reagents at 2-8 °C.

3. All reagents should equilibrate at room temperature (20-25°C) before use, and return to 2-8 °C following use.

4. Handle all materials and reagents according to the Good Laboratory Practices.

5. Serum samples, as well as positive and negative controls, should be considered potentially infectious and, therefore, the kit components that contact them should be handled with gloves and disposed as pathological waste according to legal regulations. The stop solution contains hydrochloric acid (corrosive). Skin contact generates irritation and should be handled with gloves. The rest of the kit components does not represent any potential risk for health and the environment.

6. Do not use the kit beyond expiration date or if the kit components were not kept in the conditions indicated above.

7. Do not mix components or instruction manuals from different test kits batches.

8. Care should be taken to avoid contamination of kit components.

9. If multichannel pipette containers are reused, they should be washed thoroughly with high quality water for reuse and identified to be used with the same solution.

10. Use different pipettes tips for each reagent and for each sample.

11. Include positive and negative controls, and the blank sample in duplicate in each test.

12. For reconstitution of the reagents use only distilled (Type II) or ultrapure (Type I) or high quality purified water.

Reagents preparation for one microplate (96 reactions)

a. Wash solution 1X

Prepare 400 ml of 1X Wash Solution by adding 40 ml of 10X Wash Solution to 360 ml of distilled (Type II) or ultrapure (Type I) water. Sufficient volume to process one entire microplate including the dilution of the conjugate and washes. Mix well before using. Once prepared, the 1X Washing Solution can be stored at 2-8°C for 7 days.

Note: If crystals appear in the 10X Wash solution, the bottle must be heated not above 45°C and shake well until completely dissolve.

It is recommended to prepare sufficient quantity to process the plates of the day.

b. Conjugate

Dilute 1/10 using the 1X Wash Solution. For example, to process the entire plate add 9 ml of the 1X Wash Solution and 2 ml of the 10X IgG conjugate. Homogenize very well.

Important:

-dilute the conjugate only immediately before use.

-only prepare enough volume of diluted conjugate to process the required number of wells.

-the remaining solution of the diluted conjugate should be discarded.

Preparation of controls for one microplate (96 reactions)

a. Positive and negative controls

In 1.5 ml tubes add 5 µl of each control to 495 µl of Sample diluent (1/100 dilution) and mix well. Add 25 µl of the above dilution to 75 µl of Sample diluent (dilution 1/4, final dilution factor of controls = 1/400) previously loaded into the corresponding wells and mix well (up and down with the micropipette at least 5 times). In each microplate the positive and negative controls are added at least in duplicate.

Important: shake the Sample diluent well before use.

b. Blank sample

On each plate a blank sample consisting of Sample diluent should be added. Dispense directly and in duplicate 100 µl per well of Sample diluent.

Sample preparation

Serum samples should be diluted 1/400 using the Sample diluent. In 1.5 ml tubes add 5 µl of each sample to 495 µl of Sample diluent (1/100 dilution) and mix well. Add 25 µl of the above dilution to 75 µl of Sample diluent (dilution 1/4, final dilution factor of samples = 1/400) previously loaded into the corresponding wells and mix well (up and down with the micropipette at least 5 times).

Notes:

- Fresh, chilled or frozen serum samples, free of turbidity, can be analyzed. Samples can be stored in refrigerator for 1 or 2 days. For longer periods store at -20°C and, in these cases, the samples should be completely thawed, bring to room temperature and homogenize before performing the analysis.

-Samples can be analyzed in a single well. However, it is recommended to analyze the samples at least in duplicate.

-Samples contaminated and/or not well conserved may produce erroneous results.

Instructions for the analysis of one microplate (96 reactions)

1. Allow the reagents equilibrate to room temperature (20-25°C). Do not remove the plate from its pouch until it has reached room temperature. If the entire plate is not used, separate only the 8-well strips needed to analyze the controls and samples. Store the rest of the strips, along with the desiccants, in the original foil bag (seal the open side with a couple of folds and adhesive tape) and place it in the ziplock-type polythene bag included in the kit and store in the refrigerator (2-8°C).

Note: once the original packaging of the plate has been opened, the strips can be kept in the closed polyethylene bag and in the refrigerator (2-8 °C) up to 90 days.

2. Process the controls (positive control, negative control and blank sample) as indicated in the Item Preparation of controls.

3. Process the samples as indicated in the item Sample preparation.

Note: samples can be processed in a single well or in duplicate. However, to obtain more reliable results it is recommended to run samples at least in duplicate.

4. Cover the plate and incubate at room temperature (20-25°C) with gentle orbital shaking for 1 hour.

5. Wash/rinse the plates 5 times with 200 to 300 µl of 1X Wash Solution per well. See item Reagents preparation.

Note: Remove the liquid completely after each washing by dump and tapping onto absorbent material or aspiration with pipette. Avoid plate drying between plate washings and before adding the next reagent.

6. Add 100 µl of Conjugate to each well. See item Preparation of the conjugate. Cover the plate and incubate at room temperature (20-25°C) with gentle orbital shaking for 1 hour.

7. Repeat step 5.

8. Add 100 µl of Substrate-chromogen solution to each well, and incubate with gentle orbital shaking for 10 minutes (± 1 min) at room temperature (20-25°C) and in the dark (for example, cover the plate with aluminum foil). Start counting the time after filling the first well.

9. Stop the reaction by adding 100 µl of Stop solution into each well and mix thoroughly by lightly tapping the edges of the plate. This solution must be added in the same order and speed as the Substrate-chromogen solution (see Step 8).

10. Measure the absorbance at 450 nm using a plate spectrophotometer. There is no need to perform a blank reading. Perform the reading within 15 minutes after adding the stop solution.

Note: samples with high reactivity values can lead to the appearance of a black precipitate due to the precipitation of the substrate.

Calculations

Expression of results

Results are expressed as the percentage of reactivity (PR) with respect to the positive control (PC) included in each test. The absorbance values measured at 450 nm (Abs450) of each sample are related to the Abs450 value of the PC as follows:

$$PR_{\text{SAMPLE}} = \frac{\text{Abs450 Sample}^*}{\text{Average Abs450 PC}} \times 100$$

* If the sample was run in duplicate, consider the average of the Abs450 value.

Validity criteria

To confirm the validity of the test the following criteria must be met:

- The average value of the PC (Abs450) must be higher than 1.0 (Abs450 PC > 1.0).
- The Abs450 values of the controls (PC and NC) and samples duplicates must not differ by more than 20% of the corresponding average.
- The Abs450 values of the blank sample must be less than 0.1 (Abs450 Blank sample < 0.1).
- The average Abs450 value of the NC must be less than 0,15 (Average Abs450 NC < 0.15).
- The Average Abs450 PC / Average Abs450 NC ratio must be higher than 6 (Average Abs450 PC / Average Abs450 NC > 6).

If any of these criteria is not fulfilled, the test is considered invalid. In this case, the test should be repeated strictly following the instructions detailed in the instruction manual included in the kit.

Interpretation of results

The test sample results should be interpreted as follows:

Important: respect incubation times.

Serum (dilution 1/400)

PR	Interpretation
≤ 32%	Negative
> 32% ≤ 51%	Indeterminate
> 51%	Positive

In the case of obtaining a indeterminate result, the sample should be re-tested. If the result is still indeterminate, it is recommended to test a second sample obtained after a period of at least three weeks. If the new sample is indeterminate again, the epidemiological situation should be considered and the sample should be analyzed with another technique if available.

To facilitate the calculations, the evaluation of the validity criteria of the test and the interpretation of the results, a calculation sheet is available on the company's website (<https://www.chemtest.net/elisa-caninos.php>) or request it via email to info@chemtest.net. **Important:** always check that the version number of the instruction manual provided in the kit matches the version of the spreadsheet.

Sensitivity and specificity

More than 1,300 serum samples obtained from infected animals (positive samples) and from healthy dogs (negative samples), were analyzed to determine the diagnostic sensitivity (Se) and specificity (Sp) of the assay. Based on these results, a ROC (receiver-operating-analysis) curve analysis was performed. This analysis allowed us to select the cut-off value that concurrently optimized the Se and Sp of the assay and the cut-off values for which a 100 % Se or Sp was achieved.

Table 1

Sensitivity, specificity and Youden's index.

Cut-off (PR) ^a	Se (%) ^b	Sp (%) ^b	J ^c
32	100 (95.1-100)	97.8 (95.0-99.3)	0.978
51	98.6 (92.6-99.9)	100 (98.4-100)	0.986

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Chemtest Argentina S.A.