# CHEMLIS® Brucella Glyco-iELISA

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**CHEMLIS®  
Brucella Glyco-iELISA**

## Kit components / Presentations

Components	Description	15-CLO2-1P Presentation 1 microplate	15-CLO2-2P Presentation 2 microplates	15-CLO2-4P Presentation 4 microplates
Microplate	96-well microtitre plate (12 x 8 strip wells with strip holder) coated with the antigen, sealed and stored dry	1 microplate (96 tests)	2 microplates (192 tests)	4 microplates (384 tests)
Conjugate 10X IgG	10X concentrated (horseradish peroxidase conjugated anti-bovine IgG antibodies)	1x 1.1 ml	1x 2.2 ml	1x 4.4 ml
Substrate-chromogen solution	Ready to use [3,3',5,5'-Tetramethylbenzidine in substrate buffer containing hydrogen peroxide (H <sub>2</sub> O <sub>2</sub> )]- STORE IN THE DARK	1x 15 ml	1x 30 ml	1x 60 ml
Stop solution	Ready to use (contains HCl 1%) - CORROSIVE	1x 15 ml	1x 30 ml	1x 60 ml
Wash solution	10X concentrated	1x 60 ml, 10X (sufficient for 0.6 liters)	1x 120 ml, 10X (sufficient for 1.2 liters)	2x 60 ml, 10X (sufficient for 1.2 liters)
Sample diluent	To reconstitute - CONTAINS PRESERVATIVE [Thimerosal (C <sub>9</sub> H <sub>9</sub> HgNaO <sub>2</sub> S) 0.01%].	1x 120 ml	2x 120 ml	3x 120 ml
Positive control	Inactivated serum – CONTAINS PRESERVATIVE [Thimerosal (C <sub>9</sub> H <sub>9</sub> HgNaO <sub>2</sub> S) 0.01%]	1x 0.1 ml	1x 0.1 ml	1x 0.1 ml
Negative control	Inactivated serum – CONTAINS PRESERVATIVE [Thimerosal (C <sub>9</sub> H <sub>9</sub> HgNaO <sub>2</sub> S) 0.01%]	1x 0.1 ml	1x 0.1 ml	1x 0.1 ml
Polyethylene bag	With hermetic closure. Reusable.	1 unit	1 unit	1 unit
Instruction manual		1 unit	1 unit	1 unit

## Name and application

### CHEMLIS® Brucella Glyco-iELISA

Indirect ELISA kit for the detection of antibodies against *Brucella abortus*, *Brucella melitensis* and *Brucella suis* in human serum and plasma.

CHEMLIS® Brucella Glyco-iELISA is the first indirect and recombinant ELISA kit (Enzyme Linked Immunosorbent Assay) for the detection of IgG specific antibodies against the O polysaccharide of the lipopolysaccharide (LPS) of *Brucella abortus*, *Brucella melitensis* and *Brucella suis* in serum samples or plasma obtained with heparin, sodium citrate or EDTA. Thanks to its exclusive GlycoEng technology, the kit has an excellent diagnostic performance and minimizes cross-reaction with other Gram-negative bacteria. CHEMLIS® Brucella Glyco-iELISA can be used as a screening and/or a confirmatory test.

### Introduction

Human brucellosis is one of the world's most widespread bacterial zoonoses and, over the past decade, new foci of the disease have emerged. Every year, more than 500,000 new cases are reported globally with annual incidence rates that varies widely (from < 2 to > 500 per 1,000,000 population) among different regions. The etiological agents of brucellosis are Gram-negative bacteria of the genus *Brucella*. *B. melitensis*, *B. abortus* and *B. suis* are the three main human brucellosis pathogens whose preferred natural host animals are sheep and goats, cattle, and swine, respectively. The infection is primarily transmitted by consumption of unpasteurized dairy products, direct contact with infected animals (slaughterhouse workers and veterinarians), and handling of cultures and clinical specimens in clinical, research, and production laboratories.

Brucellosis can lead to serious complications in affected patients with an important public health issue.

The wide spectrum of clinical manifestations and lack of pathognomonic symptoms make human brucellosis difficult to clinically diagnose and distinguish from several febrile conditions that often occur in the same areas; therefore, laboratory tests are essential for diagnosing the disease. Even though isolation of *Brucella*, mainly from blood, continues to be the gold standard method for the diagnosis of brucellosis, it presents several drawbacks. The slow growth of *Brucella* in *primo-cultures* delays diagnosis for several days and the sensitivity is often low ranging from 50 to 90% depending on the stage of the disease, previous use of antibiotics, the clinical specimen, and the culture method employed. Furthermore, culturing *Brucella* sp. is hazardous and not all the laboratories have suitable culture facilities. Hence, when the disease cannot be confirmed by culture, serological tests play a major role in the routine diagnosis of brucellosis.

Currently, the diagnostic methods most commonly used for human serological testing are agglutination, complement fixation, and enzyme immunoassay test. All these assays use as antigen the whole bacteria or bacterial extracts containing high concentrations of the smooth lipopolysaccharide (sLPS) since the humoral immune response to smooth brucellae is dominated by antibodies to the O-polysaccharide section of *Brucella* sLPS.

### Principles of the technique

CHEMLIS® Brucella Glyco-iELISA kit is the first indirect solid phase enzyme immunoassay based exclusively on the detection of *Brucella abortus*, *Brucella melitensis* and *Brucella suis* anti-O polysaccharide IgG antibodies. A positive result by CHEMLIS® Brucella Glyco-iELISA indicates the presence of anti-O polysaccharide antibodies. These antibodies are an excellent specific mark-

er of Brucella infection, which allows to confirm the diagnosis of infection eliminating the risk of false-positive cross-reactions due to the presence of antibodies directed against common epitopes present in the lipid A and the core of the LPS from other bacteria.

In this assay the antibodies present in the sample react with the antigen (purified recombinant glycoprotein) that coats the wells. If the sample contains antibodies against the O polysaccharide of Brucella abortus, Brucella melitensis and/or Brucella suis, they bind to the well. Upon addition of the horseradish peroxidase (HRP) conjugated antibody, a complex is formed with the IgG antibodies attached to the antigen. Free antibodies are removed by washing before adding the substrate-chromogen solution. The blue color that appears is the product of the reaction of the substrate-chromogen with the conjugate and color intensity is directly proportional to the amount of antibodies present in the sample reacting with the antigens. The higher the antibody titer present in the sample, the more intense color develops in the test wells. Adding the stop solution, upon which a color change to yellow is observed, stops the reaction. The results are measured in a microplate spectrophotometer determining the absorbance (Abs) at 450 nm.

### **Storage and expiration**

Storage: between 2 and 8°C.

Transport temperature: 4 to 15°C for 72 hours.

Proper storage of the kit and all reagents will allow the use until the expiration date.

Shelf life: 12 months. Do not use reagents from different lots or after the expiration date.

### **Materials needed but not provided**

-High precision multichannel pipette and micropipettes (10, 100 and 1000 µl).

-Containers for multichannel pipette.

-Disposable pipette tips.

-Microplate reader for absorbance determination at 450 nm.

-Orbital rotator.

-Timer.

-Tubes for the dilutions of the controls and samples.

-Container (graduated cylinder) for the preparation of the wash solution.

-Distilled water (Type II) or ultrapure water (Type I).

-Microplate covers (lid, aluminium foil or adhesive).

### **Precautions**

1. Carefully read and strictly follow all instructions detailed in the leaflet included in the kit.

2. Store the kit and all reagents at 2-8 °C.

3. All reagents should equilibrate at room temperature (20-25°C) before use, and return to 2-8 °C following use.

4. Handle all materials and reagents according to the Good Laboratory Practices.

5. The samples, as well as positive and negative controls, should be considered potentially infectious and, therefore, the kit components that contact them should be handled with gloves and disposed as pathological waste according to legal regulations. The stop solution contains hydrochloric acid (corrosive). Skin contact generates irritation and should be handled with gloves. The rest of the kit components does not represent any potential risk for health and the environment.

6. Do not use the kit beyond expiration date or if the kit components were not kept in the conditions indicated above.

7. Do not mix components or instruction manuals from different test kits batches.

8. Care should be taken to avoid contamination of kit components.

9. If multichannel pipette containers are reused, they should be washed thoroughly with high quality water for reuse

and identified to be used with the same solution.

10. Use different pipettes tips for each reagent and for each sample.

11. Include positive and negative controls, and the target sample in duplicate in each test.

12. For reconstitution of the reagents use only distilled water (Type II) or ultrapure water (Type I).

## Preparation of reagents

### a. Wash solution 1X

The 10X concentrated Wash solution should be brought to room temperature (20-25°C) and mixed to ensure dissolution of any precipitated salts.

The 10X concentrated Wash solution must be diluted 1/10 with distilled (Type II) or ultrapure (Type I) water; for example, adding 30 ml of 10X Wash solution to 270 ml of water (sufficient quantity of Wash solution 1X to process an entire microplate including the dilution of the conjugate, reconstitution of the Sample diluent and washes). Mix very well before using. Once prepared, the 1X Wash solution can be stored at 2-8 ° C for up to 7 days.

It is recommended to prepare sufficient quantity to process the plates of the day.

### b. Conjugate

Dilute 1/10 using the 1X Wash Solution. For example, to process the entire plate add 9 ml of the 1X Wash Solution and 1 ml of the 10X IgG conjugate. Homogenize very well.

#### Important:

-dilute the conjugate only immediately before use.

-only prepare enough volume of diluted conjugate to process the required number of wells.

-the remaining solution of the diluted conjugate should be discarded.

### c. Sample diluent

Reconstitute the Sample diluent by adding the 1X Wash solution previously prepared (see previous item) up to the level indicated in the corresponding bottle (120 ml) and mix very well. Prepare immediately before use.

The remaining Sample diluent can be stored in the refrigerator (2 to 8°C) up to 3 months or at -20°C for 12 months. Thaw it completely, bring to room temperature and mix very well until complete dissolution before using.

## Preparation of controls

### a. Positive (PC) and negative (NC) controls

Controls (PC and NC) should be diluted 1:400 using the Sample diluent.

Add 2 µl of the corresponding control to 800 µl of Sample diluent (dilution 1:400) and mix well. Use 100 µl of each diluted control per well and in duplicate.

### b. Sample blank (SB)

On each plate a sample blank consisting of Sample diluent should be added. Dispense directly and in duplicate 100 µl per well of Sample diluent.

## Sample preparation

Samples should be diluted 1:400 using the Sample diluent.

Add 2 µl of the sample (serum or plasma) to 800 µl of Sample diluent (dilution 1:400) and mix well. Use 100 µl of the diluted sample per well.

Notes:

-It is recommended to analyze serum or plasma (obtained using heparin, sodium citrate or EDTA as anticoagulant) samples recently obtained (fresh) and free of turbidity.

-If the samples contain particles or a high content of triglycerides, they should be clarified by centrifugation.

-If the samples are not analyzed immediately, they can be kept refrigerated (2 to 8 °C) for up to 72 hs. For longer storage, store at -20°C (-15 to -20°C), taking into account that during storage and as a consequence of the freezing and thawing process, the antibody titer may

change. If samples are frozen, they must be completely thawed, homogenized and centrifuged before testing.

-It is recommended to analyze the samples at least in duplicate.

-Samples contaminated and/or not well conserved may produce erroneous results.

## Test procedure

1. Allow the reagents equilibrate to room temperature (20-25°C). Do not remove the plate from its pouch until it has reached room temperature. If the entire plate is not used, separate only the 8-well strips needed to analyze the controls and samples. Store the rest of the strips, along with the desiccants, in the original foil bag (seal the open side with a couple of folds and adhesive tape) and place it in the ziplock-type polythene bag included in the kit and store in the refrigerator (2-8°C).

Note: once the original packaging of the plate has been opened, the strips can be stored in the hermetically sealed polyethylene bag with the desiccants, and in the refrigerator (2-8°C) for up to 4 months, as long as the indicated expiration date is not exceeded.

2. Add in duplicate 100 µl per well of the positive control, negative control and sample blank. See Item Preparation of controls.

3. Add 100 µl per well of each serum sample. See Item Sample preparation.

4. Cover the plate and incubate at room temperature (20-25°C) with gentle orbital shaking for 1 hour.

5. Wash/rinse the plates 4 times with 200 µl of 1X Wash Solution per well. See item Reagents preparation.

Note: Remove the liquid completely after each washing by dump and tapping onto absorbent material or aspiration with pipette. Avoid plate drying between plate washings and before adding the next reagent.

6. Add 100 µl of Conjugate to each well. See item Preparation of the conjugate. Cover the plate and incubate at room temperature (20-25°C) with gentle orbital shaking for 1 hour.

7. Repeat step 5.

8. Add 100 µl of Substrate-chromogen solu-

tion to each well, and incubated with gentle orbital shaking for 10 minutes ( $\pm 1$  min) at room temperature (20-25°C) and in the dark (for example, cover the plate with aluminum foil). Start counting the time after filling the first well.

9. Stop the reaction by adding 100 µl of Stop solution into each well and mix thoroughly by lightly tapping the edges of the plate. This solution must be added in the same order and speed as the Substrate-chromogen solution (see Step 8).

10. Measure the absorbance at 450 nm using a plate spectrophotometer. There is no need to perform a blank reading. Perform the reading within 15 minutes after adding the stop solution.

Note: samples with high reactivity values can lead to the appearance of a black precipitate due to the precipitation of the substrate.

**Important:** respect incubation times.

## Calculations

### Expression of results

Results are expressed as the percentage of reactivity ( $PR_{PC} = 100\%$ ) with respect to the positive control included in each test. The absorbance values measured at 450 nm ( $Abs_{450}$ ) of each sample are related to the  $Abs_{450}$  value of the positive control (PC) as follows:

$$PR_{SAMPLE} = \frac{Abs_{450} \text{ Sample}^*}{Average Abs_{450} PC} \times 100$$

\* If the sample was run in duplicate, consider the average of the  $Abs_{450}$  value.

### Validity criteria

To confirm the validity of the test the following criteria must be met:

- The average value of the PC ( $Abs_{450}$ ) must be higher than 1.6 ( $Abs_{450} PC > 1.6$ ).
- The  $Abs_{450}$  values of the controls (PC and NC) and samples duplicates must not differ by more than 20% of the corresponding average.

- The Abs<sub>450</sub> values of the blank sample must be less than 0.15 (Abs<sub>450</sub> Blank sample < 0.15).
- The average Abs<sub>450</sub> value of the NC must be less than 0.25 (Average Abs<sub>450</sub> NC < 0.25).
- The Average Abs<sub>450</sub> PC / Average Abs<sub>450</sub> NC ratio must be higher than 6.4 (Average Abs<sub>450</sub> PC / Average Abs<sub>450</sub> NC > 6.4).

If any of these criteria is not fulfilled, the test is considered invalid. In this case, the test should be repeated strictly following the instructions detailed in the Instruction manual included in the kit.

### Interpretation of results

The test sample results should be interpreted as follows:

Sample (dilution 1:400)

PR	Interpretation
≤ 35%	Non-Reactive
> 35 ≤ 45%	Indeterminate
> 45%	Reactive

In the case of obtaining a "Indeterminate" result, the sample should be re-tested. If the result is still indeterminate, it is recommended to test a second sample obtained after a period of at least three weeks. If the new sample is indeterminate again, the epidemiological situation should be considered and the sample should be analyzed with another technique if available.

To facilitate the calculations, the evaluation of the validity criteria of the test and the interpretation of the results, a calculation sheet is available on the company's website (<https://www.chemtest.net/chemlis-brucella-glyco.php>) or request it via email to [info@chemtest.net](mailto:info@chemtest.net). **Important:** always check that the version number of the instruction manual provided in the kit matches the version of the spreadsheet.

### Diagnostic sensitivity and specificity

More than 1000 serum samples obtained from patients with brucellosis caused by *B. abortus*, *B. melitensis* or *B. suis*, and from healthy donors were analyzed to determine the diagnostic sensitivity (Se) and specificity (Sp) of the assay. Based on these results, a ROC (receiver-operating-analysis) curve analysis was performed. This analysis allowed us to select the cut-off value that concurrently optimized the Se and Sp of the assay and the cut-off values for which a 100% Se or Sp was achieved. The cut-off value that simultaneously maximized the Se and Sp (maximum value of the Youden's index, J) matches the cut-off value for which the Se is 100%.

Table 1

Sensitivity, specificity and Youden's index.

Valor de corte (PR) <sup>a</sup>	Se (%) <sup>b</sup>	Sp (%) <sup>b</sup>	J <sup>c</sup>
35	100 (95.3 - 100)	98 (92.8 - 99.8)	<b>0.980</b>
45	96 (88.9 - 99.2)	100 (96.3 - 100)	0.960

(a) PR, percentage of reactivity.

(b) Se, sensitivity (TP/TP+FN) x100; Sp, specificity (TN/TN+FP) x100. Values in parentheses indicate the 95% confidence interval. TP, true positive; TN, true negative; FP, false positive; FN, false negative.

(c) J, Youden's index (Se+Sp-1).

### For technical assistance

Tel: (54-11) 5353-6066  
[info@chemtest.net](mailto:info@chemtest.net)  
[chemtest.net](http://chemtest.net)

### References

-Andrés E. Ciocchini, Diego A. Rey Serantes, Luciano J. Melli, Jeremy A. Iwashkiw, Bettina Deodato, Jorge Wallach, Mario F. Feldman, Juan E. Ugalde and Diego J. Comerchi.  
 Development and Validation of a Novel Diagnostic Test for Human Brucellosis Using a Glyco-engineered Antigen Coupled to Magnetic Beads.  
 PLoS Negl Trop Dis. 2013 Feb;7(2):e2048.  
 doi:10.1371/journal.pntd.0002048.

## Indice de símbolos / Symbol index

	Fecha de vencimiento / Expiration date
	Número de catálogo / Catalog number
	Lote / Batch code
	Para diagnóstico in vitro / For in vitro diagnostic use only
	Límites de temperatura / Temperature limits
	Consultar instrucciones para su uso / See instructions for use
	Número de determinaciones / Number of determinations
	Fabricante / Manufacturer
	Riesgo biológico / Biological risk

ELISA

Humanos

Chemtest Argentina S. A.  
Av. 25 de Mayo 1021, San Martín (C.P. 1650), Buenos Aires, Argentina.  
[info@chemtest.net](mailto:info@chemtest.net) - [chemtest.net](http://chemtest.net)