CHEMLIS® Chagas R-iELISA

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Electrochemical magnetic microbeads-based biosensor for point-of-care serodiagnosis of infectious diseases.

Biosensors & Bioelectronics; 80; 6-2016; 24-33. <https://doi.org/10.1016/j.bios.2016.01.021>.

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Tel: (54-11) 2033-1455 (Ext. 6028)

info@chemtest.net

chemtest.net

CHEMLIS®
Chagas R-iELISA

Kit components

<i>Microplate</i>	96-well microtitre plate coated with the antigens, sealed and stored dry	2 microplates (192 tests)
<i>Conjugate</i>	Lyophilized (horseradish peroxidase conjugated anti-bovine IgG antibodies)	2 units
<i>Substrate-chromogen solution</i>	Ready to use [3,3',5,5'-Tetramethylbencidine in substrate buffer containing hydrogen peroxide (H ₂ O ₂)]- STORE IN THE DARK	1 x 30 ml
<i>Stop solution</i>	Ready to use (contains HCl 1%) – CORROSIVE	1 x 30 ml
<i>Wash solution</i>	10X concentrate	1 x 120 ml, 10X (sufficient for 1.2 liters)
<i>Sample diluent</i>	To reconstitute - CONTAINS PRESERVATIVE [Thimerosal (C ₉ H ₉ HgNaO ₂ S) 0.01%].	2 x 120 ml
<i>Positive control</i>	Serum – CONTAINS PRESERVATIVE [Thimerosal (C ₉ H ₉ HgNaO ₂ S) 0.01%]	1 x 0.08 ml
<i>Negative control</i>	Serum – CONTAINS PRESERVATIVE [Thimerosal (C ₉ H ₉ HgNaO ₂ S) 0.01%]	1 x 0.08 ml
<i>Instruction manual</i>		1

Name and application

CHEMLIS® Chagas R-iELISA

Indirect ELISA kit for the detection of antibodies against *Trypanosoma cruzi* in human serum samples.

CHEMLIS® Chagas R-iELISA is an indirect ELISA kit (Enzyme Linked Immunosorbent Assay) for the detection of IgG specific antibodies against *Trypanosoma cruzi*, causative agent of Chagas disease, in human serum. Thanks to its exclusive combination of four recombinant antigens, the kit has an excellent diagnostic performance. CHEMLIS® Chagas R-iELISA can be used as a screening and/or a confirmatory test.

Introduction

American Trypanosomiasis or Chagas disease is a major health and economic problem in Latin America caused by the protozoan parasite *Trypanosoma cruzi*. This parasite is found throughout the American continent in a variety of wild and domestic mammalian reservoirs and transmitted by the bite of infected, hematophagous triatomine insects. Apart from vectorial transmission, humans can become infected with *T. cruzi* by ingestion of contaminated food and fluids, from mother to child during pregnancy and through contaminated blood transfusion or organ transplantation. It is estimated that 10-12 million people are currently infected and that up to 120 million individuals living in endemic areas are at risk of infection, and thus at risk of developing cardiac or gut pathology normally associated with chronic Chagas Disease. Increasing travel and immigration has also brought the risk of *T. cruzi* infection into non-endemic countries such as Spain, Australia and the U.S., where some autochthonous cases of transmission have already been reported. Diagnosis of Chagas disease is challeng-

ing because it is often asymptomatic during the acute phase, and evolves into a chronic stage that has different clinical forms. In addition, and due to a major decline in parasitemia during the chronic phase of the disease, *T. cruzi* detection in blood samples by direct examination, blood-culture or xenodiagnosis is difficult as well as labor- and time-consuming. Therefore, detection of anti-*T. cruzi* antibodies in body fluids (typically blood samples) is still the most effective method for demonstrating direct exposure to the parasite. At present, the most widely used serologic methods are indirect hemagglutination assay (IHA), indirect immunofluorescence (IIF), and enzyme-linked immunosorbent assay (ELISA) using total parasite homogenates or semi-purified antigenic fractions. In spite of their simplicity and low cost, these tests show variations in their reproducibility and reliability that can be attributed to the poor standardization of the antigens. The advent of recombinant DNA and peptide synthesis technologies allowed the production and one-step purification of large amounts of highly pure *T. cruzi* immunodominant antigens. An additional advantage of both recombinant antigens and synthetic peptides is that they minimize specificity problems, one of the major drawbacks of immunodiagnosis of Chagas disease.

Principles of the technique

CHEMLIS® Chagas R-iELISA kit is a solid phase enzyme immunoassay for the detection IgG antibodies against *Trypanosoma cruzi* in human serum samples. In this assay the antibodies present in the sample react with the recombinant antigens (exclusive combination of four recombinant proteins) that coat the wells. If the sample contains antibodies against these *Trypanosoma cruzi* proteins, they bind to the well. Upon addition of the horseradish peroxidase (HRP) conjugated antibody, a complex is formed with the IgG antibodies at-

tached to the antigen. Free antibodies are removed by washing before adding the substrate-chromogen solution. The blue color that appears is the product of the reaction of the substrate-chromogen with the conjugate and color intensity is directly proportional to the amount of antibodies present in the sample reacting with the antigens. The higher the antibody titer present in the sample, the more intense color develops in the test wells. Adding the stop solution, upon which a color change to yellow is observed, stops the reaction. The results are measured in a microplate spectrophotometer determining the absorbance (Abs) at 450 nm.

Kit components

-Microplate

96-well microtitre plate (12 x 8 strip wells with strip holder) coated with the antigens, sealed and stored dry.

-Wash solution

10X concentrate.

-Conjugate

Liophilized (horseradish peroxidase conjugated anti-bovine IgG antibodies).

-Substrate-chromogen solution

Ready to use [3,3',5,5'-Tetramethylbenzidine in substrate buffer containing hydrogen peroxide (H₂O₂)]- STORE IN THE DARK.

-Stop solution

Ready to use (contains HCl 1%) – CORROSIVE.

-Sample diluent

To reconstitute - CONTAINS PRESERVATIVE [Thimerosal (C₉H₉HgNaO₂S) 0.01%].

-Positive Control

Serum – CONTAINS PRESERVATIVE [Thimerosal (C₉H₉HgNaO₂S) 0.01%].

-Negative Control

Serum – CONTAINS PRESERVATIVE [Thimerosal (C₉H₉HgNaO₂S) 0.01%].

-Instruction manual

Storage and expiration

Store at room temperature between 2 and 8°C. Transport temperature: 4 to

15°C for 72 hours. Proper storage of the kit and all reagents will allow the use until the expiration date.

Expiration: 12 months. Do not use reagents from different lots or after the expiration date.

Materials needed but not provided

-High precision multichannel pipette and micropipettes (10, 100 and 1000 µl).

-Containers for multichannel pipette.

-Disposable pipette tips.

-Microplate reader for absorbance determination at 450 nm.

-Orbital rotator.

-Timer.

-Tubes for the dilutions of the controls and samples.

-1 liter container for the preparation of the wash solution.

-Distilled water (Type II) or ultrapure water (Type I).

-Microplate covers (lid, aluminium foil or adhesive).

Precautions

1. Carefully read and strictly follow all instructions detailed in the leaflet included in the kit.

2. Store the kit and all reagents at 2-8 °C.

3. All reagents should equilibrate at room temperature (20-25°C) before use, and return to 2-8 °C following use.

4. Handle all materials and reagents according to the Good Laboratory Practices.

5. Serum samples, as well as positive and negative controls, should be considered potentially infectious and, therefore, the kit components that contact them should be handled with gloves and disposed as pathological waste according to legal regulations. The stop solution contains hydrochloric acid (corrosive). Skin contact generates irritation and should be handled with gloves. The rest of the kit components does not represent any potential risk for health and the environment.

6. Do not use the kit beyond expiration date or if the kit components were not kept in the conditions indicated above.

7. Do not mix components or instruction manuals from different test kits batches.

8. Care should be taken to avoid contamination of kit components.

9. If multichannel pipette containers are reused, they should be washed thoroughly with high quality water for reuse and identified to be used with the same solution.

10. Use different pipettes tips for each reagent and for each sample.

11. Include positive and negative controls, and the target sample in duplicate in each test.

12. For reconstitution of the reagents use only distilled water (Type II) or ultrapure water (Type I).

Preparation of reagents

a. Wash solution 1X

The 10X concentrated Wash solution should be brought to room temperature (20-25°C) and mixed to ensure dissolution of any precipitated salts.

The 10X concentrated Wash solution must be diluted 1/10 with distilled (Type II) or ultrapure (Type I) water; for example, adding 40 ml of 10X Wash solution to 360 ml of water (sufficient quantity of Wash solution 1X to process one microplate: 10 ml to reconstitute the IgG conjugate, 120 ml to reconstitute the Sample diluent and the rest for the washes). Mix very well before using. Once prepared, the 1X Wash solution can be stored at 2-8 °C for up to 7 days

b. Conjugate

Reconstitute the conjugate by adding 10 ml (per flask) of 1X Wash solution previously prepared (see previous item). Add the Washing solution carefully. Allow to stand for a minute and mix well. Prepare immediately before use.

For greater accuracy it is recommended to pipet the 10 ml using two calibration marks of the pipette (for example: pipetting twice from zero to 5 ml, or once from position -1 to 9 of the pipette).

The remaining conjugate can be aliquot-ed in Eppendorf type tubes and stored at -20°C (cannot be stored in the refrigerator) for up to 30 days. During this period of time each aliquot can be thaw and re-frozen one time.

c. Sample diluent

Reconstitute the Sample diluent by adding the 1X Wash solution previously prepared (see previous item) up to the level indicated in the corresponding bottle (approximately 60 ml) and mix very well. Prepare immediately before use.

The remaining Sample diluent can be stored in the refrigerator (2 to 8°C) up to 3 months or at -20°C for a year. Thaw it completely, bring to room temperature and mix very well until complete dissolution before using.

Preparation of controls

a. Positive and negative controls

In tubes of 1.5 ml add 4 µl of each control to 800 µl of Sample diluent (1:200 dilution). Mix well before using. In each microplate positive and negative controls should be added at least in duplicate (100 µl per well).

b. Sample blank (SB)

On each plate a sample blank consisting of Sample diluent should be added. Dispense directly and in duplicate 100 µl per well of Sample diluent.

Sample preparation

Serum samples should be diluted 1:00 using the Sample diluent. Add 4 µl of serum to 800 µl of Sample diluent. Use 100 µl of diluted sample per well.

Notes:

-Fresh, chilled or frozen serum samples, free of turbidity, can be analyzed. Samples can be stored in a refrigerator for 1 or 2 days. For longer periods store at -20°C and, in these cases, the samples should be completely thawed, bringing them to room temperature, and homogenized before performing the analysis.

-It is recommended to analyze the samples at least in duplicate.

-Samples contaminated and/or not well conserved may produce erroneous results.

Instructions for the analysis

1. Allow the reagents equilibrate to room temperature (20-25°C). Do not remove the plate from its pouch until it has reached room temperature.

2. Add in duplicate 100 µl per well of the positive control, negative control and sample blank. See Item Preparation of controls.

3. Add 100 µl per well of each serum sample. See Item Sample preparation.

Note: samples can be processed in a single well or in duplicate. However, to obtain more reliable results it is recommended to run samples at least in duplicates.

4. Cover the plate and incubate at room temperature (20-25°C) with gentle orbital shaking for 1 hour.

5. Wash/rinse the plates 4 times with 200 µl of 1X Wash Solution per well. See item Reagents preparation.

Note: Remove the liquid completely after each washing by dump and tapping onto absorbent material or aspiration with pipette. Avoid plate drying between plate washings and before adding the next reagent.

6. Add 100 µl of Conjugate to each well. See item Preparation of the conjugate. Cover the plate and incubate at room temperature (20-25°C) with gentle orbital shaking for 1 hour.

7. Repeat step 5.

8. Add 100 µl of Substrate-chromogen solution to each well, and incubated with gentle orbital shaking for 10 minutes (± 1 min) at room temperature (20-25°C) and in the dark (for example, cover the plate with aluminum foil). Start counting the time after filling the first well.

9. Stop the reaction by adding 100 µl of Stop solution into each well and mix thoroughly by lightly tapping the edges of the plate. This solution must be added in the same order and speed as the Substrate-chromogen solution (see Step 8).

10. Measure the absorbance at 450 nm using a plate spectrophotometer. There is no need to perform a blank reading. Perform the reading within 15 minutes after adding the stop solution.

Note: samples with high reactivity values can lead to the appearance of a black precipitate due to the precipitation of the substrate.

Calculations

Expression of results

Results are expressed as the percentage of reactivity ($PR_{PC} = 100\%$) with respect to the positive control included in each test. The absorbance values measured at 450 nm (Abs_{450}) of each sample and the negative control (NC) are related to the Abs_{450} value of the positive control (PC) as follows:

$$PR_{CN} = \frac{\text{Average } Abs_{450} \text{ NC}}{\text{Average } Abs_{450} \text{ PC}} \times 100$$

$$PR_{\text{Sample}} = \frac{Abs_{450} \text{ Sample}^*}{\text{Average } Abs_{450} \text{ PC}} \times 100$$

* If the sample was run in duplicate, consider the average of the Abs_{450} value.

Validity criteria

To confirm the validity of the test the following criteria must be met:

- The average value of the PC (Abs_{450}) must be higher than 1.6 ($Abs_{450} \text{ PC} > 1.6$).
- The Abs_{450} values of the controls (PC and NC) and samples duplicates must not differ by more than 20% of the corresponding average.
- The Abs_{450} values of the blank sample

must be less than 0.1 (Abs_{450} Blank sample < 0.1).

- The average Abs_{450} value of the NC must be less than 0,12 (Average Abs_{450} NC < 0.12).
- The Average Abs_{450} PC / Average Abs_{450} NC ratio must be higher than 13 (Average Abs_{450} PC / Average Abs_{450} NC > 13).

If any of these criteria is not fulfilled, the test is considered invalid. In this case, the test should be repeated strictly following the instructions detailed in the Instruction manual included in the kit.

Interpretation of results

The test sample results should be interpreted as follows:

Serum (dilution 1:200)

PR	Interpretation
≤ 16%	Non-Reactive
> 16% to ≤ 26%	Indeterminate
> 26%	Reactive

In the case of obtaining an "Indeterminate" result, the sample should be re-tested. If the result is still indeterminate, it is recommended to test a second sample obtained after a period of at least three weeks. If the new sample is indeterminate again, the epidemiological situation should be considered and the sample should be analyzed with another technique if available.

Diagnostic sensitivity and specificity

More than 320 serum samples obtained from patients with acute or chronic Chagas disease and from healthy donors were analyzed to determine the diagnostic sensitivity (Se) and specificity (Sp) of the assay. Based on these results, a ROC (receiver-operating-analysis) curve analysis was performed. This analysis allowed us to select the cut-off value that concurrently optimized the Se and Sp of the assay and the cut-off values for which a

100 % Se or Sp was achieved. The cut-off value that simultaneously maximized the Se and Sp (maximum value of the Youden's index, J) matches the cut-off value for which the Se is 100% (see Table 1).

Table 1

Sensitivity, specificity and Youden's index.

Cut-off (PR) ^a	Se (%) ^b	Sp (%) ^b	J ^c
16	100 (94.6-100)	99,6 (97,9-99,9)	0.996
26	97,0 (89,6-99,6)	100 (98,6-100)	0,970

(a) PR, percentage of reactivity.

(b) Se, sensitivity (TP/TP+FN) x100; Sp, specificity (TN/TN+FP) x100. Values in parentheses indicate the 95% confidence interval. TP, true positive; TN, true negative; FP, false positive; FN, false negative.

(c) J, Youden's index (Se+Sp-1).

References

- Cortina, María E.; Melli, Luciano J.; Roberti, Mariano; Mass, Mijal I.; Longinotti, Gloria; Tropea, Salvador; Lloret, Paulina; Rey Serantes, Diego; Salomón, Francisco; Lloret, Matías; Caillava, Ana J.; Restuccia, Sabrina; Altcheh, Jaime M.; Buscaglia, Carlos A.; Malatto, Laura; Ugalde, Juan E.; Fraigi, Liliana; Moína, Carlos; Ybarra, Gabriel; Ciocchini, Andrés E. and Comerci, Diego J.
Electrochemical magnetic microbeads-based biosensor for point-of-care serodiagnosis of infectious diseases.
Biosensors & Bioelectronics; 80; 6-2016; 24-33. <https://doi.org/10.1016/j.bios.2016.01.021>.

For technical assistance

Tel: (54-11) 2033-1455 (Ext. 6028)
info@chemtest.net
chemtest.net

Elisa

Humano

Chemtest S.A.

Av. 25 de Mayo 1021, San Martín (C.P. 1650), Buenos Aires, Argentina.

info@chemtest.net - chemtest.net