CHEMLIS® E. coli Combi Glyco-iELISA

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E. coli Combi Glyco-iELISA

Kit components / Presentations

Components	Description	REF: 15-CL08-2P Presentation 2 microplates	REF: 15-CL08-4P Presentation 4 microplates
<i>Microplate</i>	96-well microtitre plate coated with the antigens, sealed and stored dry	2 microplates (9 samples in duplicate or 18 in single well)	4 microplates (18 samples in duplicate or 36 in single well)
<i>Conjugate 10X IgM</i>	10X concentrated (horseradish peroxidase conjugated anti-human IgM antibodies)	1 x 1.1 ml	1 x 2.2 ml
<i>Conjugate 10X IgG</i>	10X concentrated (horseradish peroxidase conjugated anti-human IgG antibodies)	1 x 1.1 ml	1 x 2.2 ml
<i>Substrate-chromogen solution</i>	Ready to use [3,3',5,5'-Tetramethylbenzidine in substrate buffer containing hydrogen peroxide (H ₂ O ₂)]- STORE IN THE DARK	1 x 30 ml	1 x 60 ml
<i>Stop solution</i>	Ready to use (contains HCl 1%) – CORROSIVE	1 x 30 ml	1 x 60 ml
<i>Wash solution</i>	10X concentrated	1 x 120 ml, 10X (sufficient for 1.2 liters)	1 x 120 ml, 10X (sufficient for 1.2 liters)
<i>Sample diluent</i>	To reconstitute - CONTAINS PRESERVATIVE [Thimerosal (C ₉ H ₉ HgNaO ₂ S) 0.01%].	1 x 120 ml	1 x 120 ml
<i>Positive control O157</i>	Serum – CONTAINS PRESERVATIVE [Thimerosal (C ₉ H ₉ HgNaO ₂ S) 0.01%]	1 x 0.04 ml	1 x 0.04 ml
<i>Positive control O145</i>	Serum – CONTAINS PRESERVATIVE [Thimerosal (C ₉ H ₉ HgNaO ₂ S) 0.01%]	1 x 0.04 ml	1 x 0.04 ml
<i>Positive control O121</i>	Serum – CONTAINS PRESERVATIVE [Thimerosal (C ₉ H ₉ HgNaO ₂ S) 0.01%]	1 x 0.12 ml	1 x 0.12 ml
<i>Positive control O103</i>	Serum – CONTAINS PRESERVATIVE [Thimerosal (C ₉ H ₉ HgNaO ₂ S) 0.01%]	1 x 0.04 ml	1 x 0.04 ml
<i>Negative control</i>	Serum – CONTAINS PRESERVATIVE [Thimerosal (C ₉ H ₉ HgNaO ₂ S) 0.01%]	1 x 0.20 ml	1 x 0.20 ml
<i>Instruction manual</i>		1 unit	1 unit

Name and application

CHEMLIS®

E. coli Combi Glyco-iELISA

Indirect ELISA kit for the combined and simultaneous detection of IgM and IgG antibodies against Escherichia coli O157, O145, O121 and O103 in human serum.

CHEMLIS® E. coli Combi Glyco-iELISA is an indirect ELISA kit (Enzyme Linked Immunosorbent Assay) for the detection of IgM and IgG specific antibodies against the O polysaccharide of the lipopolysaccharide (LPS) of Escherichia coli O157, O145, O121 and O103 lipopolysaccharide (LPS) in human serum. Thanks to its exclusive **GlycoEng** technology, the kit has an excellent diagnostic performance, minimizes cross-reaction with other Gram-negative bacteria and it is the first serospecific assay that allows to differentiate O157, O145, O121 and O103 E. coli infections from those caused by other E. coli serogroups.

Introduction

Shiga toxin-producing Escherichia coli (STEC) is an important foodborne pathogen associated with sporadic cases and outbreaks of diarrhea, bloody diarrhea (BD) and hemolytic uremic syndrome (HUS), a disorder characterized by microangiopathic hemolytic anemia, thrombocytopenia and acute renal failure. Diarrhea associated HUS is caused mainly by particular STEC serotypes. STEC strains are characterized by the production of Shiga toxins, Shiga 1 (Stx1) and/or Shiga 2 (Stx2), and the production of these toxins is essential in the pathogenesis of BD and HUS. E. coli O157: H7 is the most prevalent STEC serotype associated with sporadic cases and outbreaks of BD and HUS in different parts of the world; however, other STEC non-O157 serogroups such as O145, O121 and O103, among

others, can cause a similar disease.

STEC associated HUS is the most common cause of acute renal failure in children worldwide and accounts for 90% of all cases of HUS in this population. Between 10 and 15% of children with a STEC infection develop HUS. The incidence rate of post-diarrheal HUS varies according to the country and Argentina shows the highest incidence worldwide with 12-14 cases per 100,000 children less than 5 years of age per year. This incidence rate is ten times higher than in other industrialized countries with more than 400 new cases of HUS per year and a mortality rate of 2 to 5%. In Argentina, HUS is endemic and constitutes one of the main causes of acute renal failure in children. In this country, more than 70% of cases are caused by STEC O157, 12% by O145, 1.8% by O121 and approximately 2% by STEC O103. A high percentage of these children (20 - 30%) have long-term renal sequelae, which determine that this disease is a very important cause of chronic renal failure and the sixth cause of renal transplants in adolescents and the ninth cause in the entire population of transplanted patients.

The diagnosis of STEC infections is based on the isolation and characterization of STEC strains, and the detection of free fecal Shiga toxin (FFStx). However, cases of HUS without STEC isolation and negative FFSTx are frequently observed. It has been shown that STEC is shed rapidly during the diarrheal prodrome and often ceases before the onset of HUS, with the consequence that the detection of the pathogen or its toxins results extremely difficult, if not impossible. For these reasons, alternative laboratory diagnostic methods have been used such as the detection of antibodies against the LPS of STEC strains. All available serologic tests use inactivated bacteria, bacterial extracts that contain high concentrations of LPS or purified LPS as antigen. Consequently, these assays show high rates of false positive reactions due to the pres-

ence of antibodies directed against epitopes present in the lipid A and the core of LPS shared by different STEC strains as well as other bacteria.

Bacteriological methods, alone or in combination with tests for the detection of Shiga toxins in stool samples, provide evidence of STEC infection in only 20 to 50% of children with a clinical diagnosis of HUS. The combination of bacteriological and/or Shiga toxins detection techniques with the serological analysis using the CHEMLIS® E. coli Combi Glyco-iELISA kit allows to drastically increase the diagnostic efficiency in patients with STEC infections (Boletín Integrado de Vigilancia, N° 329 - SE 39 - Septiembre de 2016. Dirección Nacional de Epidemiología y Análisis de la Situación de Salud. Ministerio de Salud - Presidencia de la Nación. ISSN 2422-698X). The exclusive GlycoEng technology incorporated in the CHEMLIS® E. coli Combi Glyco-iELISA kit confers maximum specificity and sensitivity for the detection of anti-O157, O145, O121 and O103 antibodies in serum. This makes possible to differentiate the infections caused by E. coli O157, O145, O121 and O103 from those caused by other E. coli serogroups relevant as causative agents of HUS.

Principles of the technique

CHEMLIS® E. coli Combi Glyco-iELISA kit is the first indirect solid phase enzyme immunoassay based exclusively on the detection of IgM and IgG specific antibodies against the O polysaccharide of E. coli O157, O145, O121 and O103. A positive result for one of the serogroups indicates the presence of antibodies against the corresponding O polysaccharide. These antibodies are an excellent specific marker of E. coli O157, O145, O121 and/or O103 infection which can confirm the diagnosis of infection by eliminating the risk of false-positive cross-reactions due to the presence of antibodies directed against common epitopes present in the lipid A and the core of the LPS of other Gram-negative bacteria and oth-

er E. coli serogroups.

In this assay, serum samples are exposed to the wells coated with the E. coli O157, E. coli O145, E. coli O121 or E. coli O103 antigens (purified recombinant glycoproteins). If the sample contains anti-O polysaccharide antibodies, they bind to the corresponding antigen in the well. Upon addition of the horseradish peroxidase (HRP) conjugated antibody, a complex is formed with the IgM or IgG antibodies attached to the antigen in the well. Free antibodies are removed by washing before adding the substrate-chromogen solution. The blue color that appears is the product of the reaction of the substrate-chromogen with the conjugate in those samples that are positive. The color intensity is directly proportional to the amount of antibodies present in the sample reacting with the antigen. The higher the antibody titer present in the sample, the more intense color develops in the test wells. Adding the stop solution, upon which a color change to yellow is observed, stops the reaction. The results are measured in a microplate spectrophotometer determining the absorbance (Abs) at 450 nm.

Kit components

-Microplate

96-well microtitre plate coated with the antigens, sealed and stored dry. Each plate is functionalized with the four antigens distributed in different wells.

-Wash solution

10X concentrated.

-Conjugate 10X anti-IgM

10X concentrated (horseradish peroxidase conjugated anti-human IgM antibodies).

-Conjugate 10x anti-IgG

10X concentrated (horseradish peroxidase conjugated anti-human IgG antibodies).

-Substrate-chromogen solution

Ready to use [3,3',5,5'-Tetramethylbenzidine in substrate buffer containing hy-

drogen peroxide (H₂O₂)]- STORE IN THE DARK.

-Stop solution

Ready to use (contains HCl 1%) – CORROSIVE.

-Sample diluent

To reconstitute - CONTAINS PRESERVATIVE [Thimerosal (C₉H₉HgNaO₂S) 0.01%].

-Positive Control O157

Serum – CONTAINS PRESERVATIVE [Thimerosal (C₉H₉HgNaO₂S) 0.01%].

-Positive Control O145

Serum – CONTAINS PRESERVATIVE [Thimerosal (C₉H₉HgNaO₂S) 0.01%].

-Positive Control O121

Serum – CONTAINS PRESERVATIVE [Thimerosal (C₉H₉HgNaO₂S) 0.01%].

-Positive Control O103

Serum – CONTAINS PRESERVATIVE [Thimerosal (C₉H₉HgNaO₂S) 0.01%].

-Negative Control

Serum – CONTAINS PRESERVATIVE [Thimerosal (C₉H₉HgNaO₂S) 0.01%].

-Instruction manual

Storage and expiration

Store at room temperature between 2 and 8°C. Transport temperature: 4 to 15°C for 72 hours. Proper storage of the kit and all reagents will allow the use until the expiration date.

Expiration: 12 months. Do not use reagents from different lots or after the expiration date.

Materials needed but not provided

-High precision multichannel pipette and micropipettes (10, 100 and 1000 µl).

-Containers for multichannel pipette.

-Disposable pipette tips.

-Microplate reader for absorbance determination at 450 nm.

-Orbital rotator.

-Timer.

-Tubes for the dilutions of the controls and samples.

-1 liter container for the preparation of the wash solution.

-Distilled water (Type II) or ultrapure water (Type I).

-Microplate covers (aluminium foil).

Precautions

1. Carefully read and strictly follow all instructions detailed in the leaflet included in the kit.

2. Store the kit and all reagents at 2-8 °C.

3. All reagents should equilibrate at room temperature (20-25°C) before use, and return to 2-8 °C following use.

4. Handle all materials and reagents according to the Good Laboratory Practices.

5. Serum samples, as well as positive and negative controls, should be considered potentially infectious and, therefore, the kit components that contact them should be handled with gloves and disposed as pathological waste according to legal regulations. The stop solution contains hydrochloric acid (corrosive). Skin contact generates irritation and should be handled with gloves. The rest of the kit components does not represent any potential risk for health and the environment.

6. Do not use the kit beyond expiration date or if the kit components were not kept in the conditions indicated above.

7. Do not mix components or instruction manuals from different test kits batches.

8. Care should be taken to avoid contamination of kit components.

9. If multichannel pipette containers are reused, they should be washed thoroughly with high quality water for reuse and identified to be used with the same solution.

10. Use different pipettes tips for each reagent and for each sample.

11. Include positive and negative controls, and the target sample in duplicate in each test.

12. For reconstitution of the reagents use only distilled water (Type II) or ultrapure water (Type I).

Preparation of reagents

a. Wash solution 1X

The 10X concentrated Wash solution should be brought to room temperature (20-25°C) and mixed to ensure dissolution of any precipitated salts.

The 10X concentrated Wash solution must be diluted 1/10 with distilled (Type II) or ultrapure (Type I) water; for example, adding 50 ml of 10X Wash solution to 450 ml of water (sufficient quantity of Wash solution 1X to process two microplates: 18 ml to dilute the IgM and IgG conjugates, 120 ml to reconstitute the Sample diluent and rest for the washes). Mix very well before using. Once prepared, the 1X Wash solution can be stored at 2-8°C for up to 7 days.

It is recommended to prepare sufficient quantity to process the plates of the day.

b. Conjugates

Dilute each conjugate 1/10 using the 1X Wash Solution. For example, to process the entire plate add 9 ml of the 1X Wash Solution and 1 ml of the 10X IgM conjugate. Proceed in the same way for the 10X IgG conjugate. Homogenize very well.

Important:

-dilute the conjugates only immediately before use.

-only prepare enough volume of diluted conjugate to process the required number of wells.

-the remaining solution of the diluted conjugate should be discarded.

c. Sample diluent

Reconstitute the Sample diluent by adding the 1X Wash solution previously prepared (see previous item) up to the level indicated in the corresponding bottle (approximately 120 ml) and mix very well. Prepare immediately before use.

The remaining Sample diluent can be stored in the refrigerator (2 to 8°C) up to 3 months or at -20°C for a year. In the latter case, thaw it completely, bring to

room temperature and mix very well until complete dissolution before using.

Preparation of controls

For one microplate according to Option 1 or for two microplates according to Option 2 (see Figure 1).

a. Positive controls (PC)

Positive controls for O157 and O103. The same dilution of the PC is used for the determination of IgM and IgG. In 1.5 ml tubes add 5 µl of each PC to 500 µl of Sample Diluent (1:100 dilution).

Positive control for O145. It is diluted 1:100 for the determination of IgM and 1:1000 for IgG.

-PC_{O145} for IgM: add 5 µl de PC_{O145} to 500 µl of Sample diluent (1:100 dilution). Mix very well.

-PC_{O145} for IgG: add 50 µl of the 1:100 dilution of the PC_{O145} previously prepared for IgM to 450 µl of the Sample diluent (final dilution: 1:1,000).

Positive control for O121. The same dilution of the PC is used for the determination of IgM and IgG. In 1.5 ml tubes add 25 µl of the PC0121 to 475 µl of Sample Diluent (1/20 dilution).

Mix very well before using. In each microplate all the positive controls are assayed in duplicate (100 µl per well) according to the scheme of Figure 1.

b. Negative control (NC)

In tubes of 2 ml (or more volume) add 20 µl of the NC to 2,000 µl of Sample diluent (1:100 dilution). Mix well before using. In each microplate the NC is assayed in duplicate (100 µl per well) for IgM and IgG, and for each antigen according to the scheme of Figure 1.

c. Sample blank (SB)

On each plate a sample blank consisting of Sample diluent should be added. Dispense directly and in duplicate 100 µl per well of Sample diluent according to the scheme of Figure 1.

Sample preparation

For one microplate according to Option 1 or for two microplates according to Option 2 (see Figure 1).

Serogroups O157, O145 and O121: serum samples should be diluted 1:100 using the Sample diluent. For example, add 10 µl of serum to 1,000 µl of Sample diluent (for single analysis) or 20 µl in 2,000 µl (for analysis in duplicate). Use 100 µl of diluted sample per well.

Serogroup O103: serum samples should be diluted 1:400 using the Sample diluent. For example, add 2.5 µl of serum to 1,000 µl of Sample diluent (for single analysis) or 5 µl in 2,000 µl (for analysis in duplicate). Use 100 µl of diluted sample per well.

Notes:

-Fresh, chilled or frozen serum samples, free of turbidity, can be analyzed. Samples can be stored in a refrigerator for 1 or 2 days. For longer periods store at -20°C and, in these cases, the samples should be completely thawed, bringing them to room temperature, and homogenized before performing the analysis.

-It is recommended to analyze the samples at least in duplicate.

-Samples contaminated and/or not well conserved may produce erroneous results.

Instructions for the analysis

1. Allow the reagents equilibrate to room temperature (20-25°C). Do not remove the plate from its pouch until it has reached room temperature.

2. Add in duplicate 100 µl per well of the positive control, negative control and sample blank. See Item Preparation of controls.

3. Add 100 µl per well of each serum sample. See Item Sample preparation.

Note: samples can be processed in a single well or in duplicate. However, to obtain more reliable results it is recommended to run samples at least in duplicates.

4. Cover the plate and incubate at room temperature (20-25°C) with gentle orbital shaking for 1 hour.

5. Wash/rinse the plates 4 times with 200 µl

of 1X Wash Solution per well. See item Reagents preparation.

Note: Remove the liquid completely after each washing by dump and tapping onto absorbent material. Avoid plate drying between plate washings and before adding the next reagent.

6. Add 100 µl of Conjugate to each well. See item Preparation of the conjugate. Cover the plate and incubate at room temperature (20-25°C) with gentle orbital shaking for 1 hour.

7. Repeat step 5.

8. Add 100 µl of Substrate-chromogen solution to each well, and incubated with gentle orbital shaking for 10 minutes (± 1 min) at room temperature (20-25°C) and in the dark (for example, cover the plate with aluminum foil). Start counting the time after filling the first well.

9. Stop the reaction by adding 100 µl of Stop solution into each well and mix thoroughly by lightly tapping the edges of the plate. This solution must be added in the same order and speed as the Substrate-chromogen solution (see Step 8).

10. Measure the absorbance at 450 nm using a plate spectrophotometer. There is no need to perform a blank reading. Perform the reading within 5 and 15 minutes after adding the stop solution.

Note: samples with high reactivity values can lead to the appearance of a black precipitate due to the precipitation of the substrate.

Important: respect incubation times.

Calculations

Expression of results

Results are expressed as the percentage of reactivity ($PR_{PC} = 100\%$) with respect to the positive control included in each test. The absorbance values measured at 450 nm (Abs_{450}) of each sample are related to the Abs_{450} value of the positive control (PC) as follows:

Determination of IgM antibodies	
	Average $Abs_{450} S_n \cdot X_{IgM}$
$PR (S_n \cdot X_{IgM}) =$	----- x 100
	Average $Abs_{450} PC \cdot X_{IgM}$

Determination of IgG antibodies	
	Average $Abs_{450} S_n \cdot X_{IgG}$
$PR (S_n \cdot X_{IgG}) =$	----- x 100
	Average $Abs_{450} PC \cdot X_{IgG}$

S, sample.

n, sample number.

X, serogroup: O157, O145, O121 or O103.

Example:

The following table shows an example of calculation for serogroup O157 with possible results for the determination of IgM antibodies.

	Abs_{450} IgM O157	PR IgM O157
Sample 1	1.20	= $(1.2/1.15) \times 100$ = 104.3%
Sample 2	0.11	= $(0.11/1.15) \times 100$ = 9.6%
PC O157	1.15	
NC O157	0.06	
SB	0.05	

Validity criteria

To confirm the validity of the IgM test the following criteria must be met:

Serogroup	Abs_{450} PC (*)	Abs_{450} NC	Abs_{450} SB	Abs_{450} PC / Abs_{450} NC
O157	> 0.9	< 0.1	< 0.1	> 9
O145	> 0.6	< 0.1	< 0.1	> 6
O121	> 0.8	< 0.1	< 0.1	> 8
O103	> 0.6	< 0.1	< 0.1	> 6

(*) In all cases consider the average value of Abs.

To confirm the validity of the IgG test the following criteria must be met:

Serogrupo	Abs_{450} CP (*)	Abs_{450} CN	Abs_{450} BM	Abs_{450} CP / Abs_{450} CN
O157	> 1.3	< 0.35	< 0.1	> 3,7
O145	> 1.0	< 0.1	< 0.1	> 10
O121	> 1.8	< 0.2	< 0.1	> 9
O103	> 1.1	< 0.7	< 0.1	> 1.6

(*) In all cases consider the average value of Abs.

-In all cases, the Abs_{450} values of the controls (PC and NC) and samples duplicate must not differ by more than 20% of the corresponding average.

-If any of these criteria is not fulfilled, the test is considered invalid. In this case, the test should be repeated strictly following the instructions detailed in the instruction manual included in the kit.

Interpretation of results

The test sample results should be interpreted as indicated in Table 1.

In the case of obtaining an "Indeterminate" result or near the cut-off value (for O121 IgG), the sample should be re-tested. If the result is still indeterminate, it is recommended to test one or more samples obtained in the following days. If with the new samples the results are still indeterminate, the clinical and epidemiological situation should be considered and the sample should be analyzed with other technique/s if available.

Diagnostic sensitivity and specificity

More than 150 serum samples obtained from children under 8 years of age were

analyzed to determine the diagnostic sensitivity (Se) and specificity (Sp) of the assay. The group of positive samples included samples obtained from patients with clinical diagnosis of Bloody Diarrhea (BD) or Hemolytic Uremic Syndrome (HUS) and a stool culture positive for STEC O157, STEC O145, STEC O121 or STEC O103. The group of negative samples included serum samples obtained from healthy children or from patients with other non-related infectious diseases.

Based on these results, a ROC (receiver-operating-analysis) curve analysis was performed. This analysis allowed us to select the cut-off value that concurrently optimized the Se and Sp of the assay and the cut-off values for which a 100 % Se or Sp was achieved (see Table 2).

For O103 the cut-off values were calculated based on the mean and the standard deviation of the reactivity values of the negative samples (samples obtained from healthy children under 12 years of age or with unrelated pathologies) (see Table 2).

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
Tel: (54-11) 2033-1455 (Ext. 6028)
info@chemtest.net
chemtest.net

Option 1: analysis of one microplate.

Option 2: analysis of two microplates.


Figure 1. Scheme of the ELISA plates for the determination of IgM and IgG antibodies. Option 1: proposed scheme for the analysis of one microplate (allows to analyze a maximum of 3 samples in duplicate or 6 samples in single well). Option 2: proposed scheme for the analysis of two microplates (allows to analyze a maximum of 9 samples in duplicate or 18 samples in single well). The positions on the microplate for the positive control of each serogroup (PC₀₁₅₇, PC₀₁₄₅, PC₀₁₂₁ and PC₀₁₀₃), negative control (NC) and sample blank (SB) are indicated. The rest of the wells are available for the analysis of the serum samples (S).

Table 1. Interpretation of the results.

Serogroup O157

Determination of IgM antibodies Serum (dilution 1:100)	
PR	Interpretation
≤ 24%	Non Reactive
> 24% to < 43%	Indeterminate
≥ 43%	Reactive

Determination of IgG antibodies Serum (dilution 1:100)	
PR	Interpretation
≤ 47%	Non Reactive
> 47% to < 60%	Indeterminate
≥ 60%	Reactive

Serogroup O145

Determination of IgM antibodies Serum (dilution 1:100)	
PR	Interpretation
≤ 18%	Non Reactive
> 18% to < 46%	Indeterminate
≥ 46%	Reactive

Determination of IgG antibodies Serum (dilution 1:100)	
PR	Interpretation
≤ 23%	Non Reactive
> 23% to < 41%	Indeterminate
≥ 41%	Reactive

Serogroup O121

Determination of IgM antibodies Serum (dilution 1:100)	
PR	Interpretation
≤ 25%	Non Reactive
> 25% to < 40%	Indeterminate
≥ 40%	Reactive

Determination of IgG antibodies Serum (dilution 1:100)	
PR	Interpretation
≤ 44%	Non Reactive
> 44%	Reactive

Serogroup O103

Determination of IgM antibodies Serum (dilution 1:400)	
PR	Interpretation
≤ 16%	Non Reactive
> 16% to < 40%	Indeterminate
≥ 40%	Reactive

Determination of IgG antibodies Serum (dilution 1:400)	
PR	Interpretation
≤ 50%	Non Reactive
> 50% to < 80%	Indeterminate
≥ 80%	Reactive

Table 2. Sensitivity, specificity and Youden's index of the assays for the selected cut-off values.

Serogroup O157				
Glyco-iELISA	Cut-off value (PR) ^a	Se (%) ^b	Sp (%) ^b	J ^c
IgM				
	> 24	100 (91.0 - 100)	94.6 (86.7 - 98.5)	0.946
	> 43	92.3 (79.1 - 98.4)	100 (95.1 - 100)	0.923
IgG				
	> 47	100 (91.0 - 100)	95.9 (88.6 - 99.2)	0.956
	> 60	94.9 (82.7 - 99.4)	100 (95.1 - 100)	0.949
Serogroup O145				
Glyco-iELISA	Cut-off value (PR) ^a	Se (%) ^b	Sp (%) ^b	J ^c
IgM				
	> 18	100 (87.2 - 100)	98.6 (92.7 - 100)	0.986
	> 46	85.2 (66.3 - 95.8)	100 (95.1 - 100)	0.852
IgG				
	> 23	100 (87.2 - 100)	82.4 (71.8 - 90.3)	0.824
	> 41	96.3 (81.0 - 99.9)	100 (95.1 - 100)	0.963
Serogroup O121				
Glyco-iELISA	Cut-off value (PR) ^a	Se (%) ^b	Sp (%) ^b	J ^c
IgM				
	> 25	100 (66.4 - 100)	97.3 (90.6 - 99.7)	0.973
	> 40	88.9 (51.8 - 99.7)	100 (95.1 - 100)	0.889
IgG				
	> 44	100 (66.4 - 100)	100 (95.1 - 100)	1.000
Serogroup O103				
Glyco-iELISA	Cut-off value (PR) ^a	Se (%) ^b	Sp (%) ^b	J ^c
IgM				
	> 16	100 (87.2 - 100)	96.3 (87.2 - 99.5)	0.963
	> 40	66.7 (9.4 - 99.2)	100 (93.3 - 100)	0.667
IgG				
	> 50	100 (29.2 - 100)	96.3 (87.2 - 99.5)	0.963
	> 80	66.7 (9.3 - 99.3)	100 (93.4 - 100)	0.667

(a) PR, percentage of reactivity.

(b) Se, sensitivity (TP/TP+FN) x100; Sp, specificity (TN/TN+FP) x100. Values in parentheses indicate the 95% confidence interval. TP, true positive; TN, true negative; FP, false positive; FN, false negative.

(c) J, Youden's index (Se+Sp-1).

ELISA

Humanos

Chemtest Argentina S.A.
Av. 25 de Mayo 1021, San Martín (C.P. 1650), Buenos Aires, Argentina.
info@chemtest.net - chemtest.net