VETLIS® Brucella Glyco-iELISA Caprinos

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-Ciocchini, AE; Rey Serantes, DA; Melli, LJ; Guidolin, LS; Iwashkiw, J; Elena, S; Franco, C; Nicola, AM; Feldman, MF; Comerci, DJ and Ugalde, JE. A bacterial engineered glycoprotein as a novel antigen for diagnosis of bovine brucellosis. Veterinary Microbiology, Aug. 2014, 27; 172(3-4):455-65.

-Andrés E. Ciocchini, Diego A. Rey Serantes, Luciano J. Melli, Jeremy A. Iwashkiw, Bettina Deodato, Jorge Wallach, Mario F. Feldman, Juan E. Ugalde and Diego J. Comerci. Development and Validation of a Novel Diagnostic Test for Human Brucellosis Using a Glyco-engineered Antigen Coupled to Magnetic Beads. *PLoS Negl Trop Dis*. 2013 Feb;7(2):e2048.

-Iwashkiw, J., Fentabil, M., Faridmoayer, A., Mills, D.C., Peppler, M. Czibener, C., Ciocchini, A.E., Comerci, D.J., Ugalde, J.E. and Feldman, M.F. Exploiting the *Campylobacter jejuni* protein glycosylation system for glycoengineering vaccines and diagnostic tools directed against brucellosis. *Microbial Cell Factories* 2012 Jan 25;11-13.

VETLIS®
Brucella Glyco-iELISA Caprine (serum)

Kit components

<i>Microplate</i>	96-well microtitre plate coated with the antigen, sealed and stored dry	2 microplates (192 tests)
<i>Conjugate</i>	Lyophilized (horseradish peroxidase conjugated anti-caprine IgG antibodies)	2 units
<i>Substrate-chromogen solution</i>	Ready to use [3,3',5,5'-Tetramethylbencidine in substrate buffer containing hydrogen peroxide (H ₂ O ₂)]- STORE IN THE DARK	1 x 30 ml
<i>Stop solution</i>	Ready to use (contains HCl 1%) – CORROSIVE	1 x 30 ml
<i>Wash solution</i>	10X concentrate	1 x 120 ml, 10X (sufficient for 1.2 liters)
<i>Caprine Positive Control Serum</i>	Serum – CONTAINS PRESERVATIVE [Thimerosal (C ₉ H ₉ HgNaO ₂ S) 0.01%]	1 x 0.025 ml
<i>Caprine Negative Control Serum</i>	Serum – CONTAINS PRESERVATIVE [Thimerosal (C ₉ H ₉ HgNaO ₂ S) 0.01%]	1 x 0.025 ml
<i>Instruction manual</i>		1

Name and application

VETLIS®

Brucella Glyco-iELISA Caprine (serum)

Indirect ELISA kit for the detection of anti *Brucella melitensis* antibodies in caprine serum.

The kit **VETLIS® Brucella Glyco-iELISA Caprine** is an ELISA kit (Enzyme Linked Immunosorbent Assay) for the indirect detection of specific antibodies against the O polysaccharide of the lipopolysaccharide (LPS) of *Brucella melitensis* in caprine serum. The kit has an excellent diagnostic performance and minimizes cross-reaction with other Gram-negative bacteria.

Introduction

Brucellosis is a highly contagious zoonosis disease caused by Gram-negative bacteria of the genus *Brucella* that affects animals and human beings. Animal brucellosis has a major economic impact because the infection causes abortions, stillbirths and reduces fertility in herds, while brucellosis in humans is a debilitating disease characterized by fever, sweating and pain that can progress to chronicity with serious health complications.

Brucella melitensis is the etiological agent of caprine brucellosis and one of the main human brucellosis pathogens together with *B. suis* and *B. abortus*. Due to the high economic impact of the disease in animal production and health, and to the risk of transmission to humans, most countries have implemented programs to control and/or eradicate brucellosis in goats. Brucellosis in humans can be debilitating and disabling and is an important public health problem. In the absence of a vaccine against human brucellosis, disease prevention depends primarily on the control of brucellosis in animals, the natural reservoir of the disease.

Control of the disease in goats depends on vaccination with the strain *B. melitensis* Rev.1, in combination with the detection and slaughter of infected animals. The gold standard method for confirmation of the infection is the isolation of the pathogen; however, the slow growth of brucellae in primary cultures (up to 7 days), the risk involved in their handling and the poor sensitivity makes diagnosis based solely on isolation of brucellae non-effective, not always feasible and expensive. For these reasons, the **laboratory diagnosis of caprine brucellosis is mainly based on the detection of specific antibodies against the infectious agent in serum samples.**

Principle of the technique

The **VETLIS® Brucella Glyco-iELISA Caprine** kit is the first indirect solid phase enzyme immunoassay based exclusively on the detection of *Brucella* anti-O polysaccharide IgG antibodies. A positive result by *Brucella Glyco-iELISA Caprine* indicates the presence of anti-O polysaccharide antibodies. These antibodies are an excellent specific marker of *brucella* infection, which can confirm the diagnosis of the infection by eliminating the risk of false-positive cross-reactions due to the presence of antibodies against common epitopes present in the lipid A and the core of the LPS from other bacteria. In this assay, serum samples are exposed to the wells coated with the *Brucella* antigen (purified recombinant glycoprotein). If the sample contains anti *Brucella melitensis* O polysaccharide antibodies, they bind to the antigen in the well. Upon addition of the horseradish peroxidase (HRP) conjugated antibody, a complex is formed with the bovine IgG antibodies attached to the antigen in the well. Free antibodies are removed by washing before adding the substrate-chromogen solution. The blue color that appears is the product of the reaction of the substrate-chromo-

gen with the conjugate in those samples that are positive. The color intensity is directly proportional to the amount of antibodies present in the sample reacting with the antigen. The higher the antibody titer present in the sample, the more intense color develops in the test wells. The reaction is stopped by adding the stop solution, upon which a color change to yellow is observed. The reading of the results is performed in a microplate spectrophotometer measuring the absorbance (Abs) at 450 nm.

Kit components

-Microplate

96-well microtitre plate coated with the antigen, sealed and stored dry.

-Wash solution

10X concentrate.

-Conjugate

Liophilized (horseradish peroxidase conjugated anti-caprine IgG antibodies).

-Substrate-chromogen solution

Ready to use [3,3',5,5'-Tetramethylbenzidine in substrate buffer containing hydrogen peroxide (H₂O₂)]- STORE IN THE DARK.

-Stop solution

Ready to use (contains HCl 1%) – CORROSIVE.

-Caprine Positive Control Serum

Serum – CONTAINS PRESERVATIVE [Thimerosal (C₉H₉HgNaO₂S) 0.01%].

-Caprine Negative Control Serum

Serum – CONTAINS PRESERVATIVE [Thimerosal (C₉H₉HgNaO₂S) 0.01%].

-Instruction manual

Storage and expiration

Store at room temperature between 2 and 8°C. Transport temperature: 4 to 15°C for 72 hours. Proper storage of the kit and all reagents will allow the use until the expiration date.

Expiration: 12 months. Do not use reagents from different lots or after the expiration date.

Materials needed but not provided

-High precision multichannel pipette and micropipettes (10, 100 and 1000 µl).

-Containers for multichannel pipette.

-Disposable pipette tips.

-Microplate reader for absorbance determination at 450 nm.

-Orbital rotator.

-Timer.

-Tubes for the dilutions of the controls and samples.

-1 liter container for the preparation of the wash solution.

-Distilled (Type II) or ultrapure (Type I) water or any similar high quality water.

-Microplate covers (lid, aluminium foil or adhesive).

Precautions

1. Carefully read and strictly follow all instructions detailed in the Instruction manual included in the kit.

2. Store the kit and all reagents at 2-8 °C.

3. All reagents should equilibrate at room temperature (20-25°C) before use, and return to 2-8 °C following use.

4. Handle all materials and reagents according to the Good Laboratory Practices.

5. Serum samples, as well as positive and negative controls, should be considered potentially infectious and, therefore, the kit components that contact them should be handled with gloves and disposed as pathological waste according to legal regulations. The stop solution contains hydrochloric acid (corrosive). Skin contact generates irritation and should be handled with gloves. The rest of the kit components does not represent any potential risk for health and the environment.

6. Do not use the kit beyond expiration date or if the kit components were not kept in the conditions indicated above.

7. Do not mix components or instruction manuals from different test kits batches.

8. Care should be taken to avoid contamination of kit components.

9. If multichannel pipette containers

are reused, they should be washed thoroughly with high quality water for reuse and identified to be used with the same solution.

10. Use different pipettes tips for each reagent and for each sample.

11. Include positive and negative controls, and the target sample in duplicate in each test.

12. For reconstitution of the reagents use only distilled (Type II) or ultrapure (Type I) water.

Reagents preparation for one microplate (96 reactions)

a. Wash solution 1X

Prepare 400 ml of 1X Wash Solution (enough for a microplate) by adding 40 ml of 10X Wash Solution to 360 ml of distilled (Type II) or ultrapure (Type I) water. Mix well before using. Once prepared, the 1X Washing Solution can be stored at 2-8°C for 7 days.

Note: If crystals appear in the 10X Wash solution, the bottle must be heated not above 45°C and shake well until completely dissolve.

b. Conjugate

Reconstitute the conjugate by adding 10 ml (per flask) of 1X Wash Solution previously prepared (see previous item). Add the Washing solution carefully. Allow to stand for a minute and mix well (cover the bottle and homogenize twice by head-tail shaking). Prepare before use. The reconstituted conjugate should be stored at -20°C (cannot be stored in the refrigerator) for not more than 14 days and during this period can be thaw and re-freeze just once.

Preparation of controls for one microplate (96 reactions)

a. Positive and negative controls

In tubes of 1.5 ml add 5 µl of each control to 1000 µl of 1X Wash Solution (1:200 dilution). Mix well before using. In each microplate positive and negative controls should be added at least in duplicate

(100 µl per well).

b. Blank sample

On each plate a blank sample consisting of 1X Wash Solution should be added. Dispense directly and in duplicate 100 µl per well of 1X Wash Solution.

Sample preparation

Serum samples should be diluted 1:200 using the 1X Wash Solution. Add 5 µl of serum to 1000 µl of 1X Wash Solution. Use 100 µl of diluted sample per well.

Notes:

-Fresh, chilled or frozen serum samples, free of turbidity, can be analyzed. Samples can be stored in refrigerator for 1 or 2 days. For longer periods store at -20°C and, in these cases, the samples should be completely thawed, bring to room temperature and homogenize before performing the analysis.

-Samples can be analyzed in a single well. However, it is recommended to analyze the samples at least in duplicate.

-Samples contaminated and/or not well conserved may produce erroneous results.

Instructions for the analysis of one microplate (96 reactions)

1. Allow the reagents equilibrate to room temperature (20-25°C). Do not remove the plate from its pouch until it has reached room temperature.

2. Add in duplicate 100 µl per well of the positive control, negative control and the blank sample. See Item Preparation of controls.

3. Add 100 µl per well of each serum sample. See the item Sample preparation.

Note: samples can be processed in a single well or in duplicate. However, to obtain more reliable results it is recommended to run samples at least in duplicates.

4. Cover the plate and incubate at room temperature (20-25°C) with gentle orbital shaking for 1 hour.

5. Wash/rinse the plates 4 times with 200 to 300 µl of 1X Wash Solution per well. See item Reagents preparation.

Note: Remove the liquid completely after each washing by dump and tapping onto absorbent

material or aspiration with pipette. Avoid plate drying between plate washings and before adding the next reagent.

6. Add 100 µl of Conjugate to each well. See item Preparation of the conjugate. Cover the plate and incubate at room temperature (20-25°C) with gentle orbital shaking for 1 hour.

7. Repeat step 5.

8. Add 100 µl of Substrate-chromogen solution to each well, and incubated with gentle orbital shaking for 10 minutes (± 1 min) at room temperature (20-25°C) and in the dark (for example, cover the plate with aluminum foil). Start counting the time after filling the first well.

9. Stop the reaction by adding 100 µl of Stop solution into each well and mix thoroughly by lightly tapping the edges of the plate. This solution must be added in the same order and speed as the Substrate-chromogen solution (see Step 8).

10. Measure the absorbance at 450 nm using a plate spectrophotometer. There is no need to perform a blank reading. Perform the reading within 15 minutes after adding the stop solution.

Note: samples with high reactivity values can lead to the appearance of a black precipitate due to the precipitation of the substrate.

Calculations

Expression of results

Results are expressed as the percentage of reactivity (PR) with respect to the positive control included in each test. The absorbance values measured at 450 nm (Abs450) of each sample and the negative control (NC) are related to the Abs450 value of the positive control (PC) as follows:

$$PR_{NC} = \frac{\text{Average Abs}_{450} \text{ NC}}{\text{Average Abs}_{450} \text{ PC}} \times 100$$

$$PR_{\text{SAMPLE}} = \frac{\text{Abs450 Sample}^*}{\text{Average Abs}_{450} \text{ PC}} \times 100$$

* If the sample was run in duplicate, consider the average of the Abs450 value.

Validity criteria

To confirm the validity of the test the following criteria must be met:

- The average value of the PC (Abs450) must be higher than 1.2 (Abs450 PC > 1.2).
- The Abs450 values of the controls (PC and NC) and samples duplicates must not differ by more than 20% of the corresponding average.
- The Abs450 values of the blank sample must be less than 0.1 (Abs450 Blank sample < 0.1).
- The average Abs450 value of the NC must be less than 0.15 (Average Abs450 NC < 0.15).
- The Average Abs450 PC / Average Abs450 NC ratio must be higher than 8 (Average Abs450 PC / Average Abs450 NC > 8).

If any of these criteria is not fulfilled, the test is considered invalid. In this case, the test should be repeated strictly following the instructions detailed in the Instruction manual included in the kit.

Interpretation of results

The test sample results should be interpreted as follows:

Serum (dilution 1:200)

PR	Interpretation
≤ 29%	Negative
> 29% a < 43%	Doubtful
≥ 43%	Positive

In the case of obtaining a doubtful result, the sample should be re-tested. If the result is still doubtful, it is recommended to test a second sample obtained after a period of at least three weeks. If the new sample is doubtful again, the epidemiological situation should be considered and the sample should be analyzed with

another technique if available.

Sensitivity and specificity

More than 400 serum samples obtained from infected animals (positive samples) as well as from animals coming from establishments free of brucellosis (negative samples) were analyzed to determine the diagnostic sensitivity (Se) and specificity (Sp) of the assay.

Based on these results, a ROC (receiver-operating-analysis) curve analysis was performed. This analysis allowed us to select the cut-off value that concurrently optimized the Se and Sp of the assay (PR=29), that coincides with the cut-off value for which the Se is 100%, and the cut-off value for which the Sp is 100% (PR=43) (Table 1).

Table 1

Sensitivity, specificity and Youden's index.

Cut-off (PR) ^a	Se (%) ^b	Sp (%) ^b	J ^c
29	100 (94.3-100)	96.4 (93.7-98.2)	0.964
43	95.2 (86.7-99.1)	100 (98.8-100)	0.952

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-Ciochini, AE; Rey Serantes, DA; Melli, LJ; Guidolin, LS; Iwashkiw, J; Elena, S; Franco, C; Nicola, AM; Feldman, MF; Comerci, DJ and Ugalde, JE. A bacterial engineered glycoprotein as a novel antigen for diagnosis of bovine brucellosis. Veterinary Microbiology, Aug. 2014, 27; 172(3-4):455-65.

-Andrés E. Ciochini, Diego A. Rey Serantes, Luciano J. Melli, Jeremy A. Iwashkiw, Bettina Deodato, Jorge Wallach, Mario F. Feldman, Juan E. Ugalde and Diego J. Comerci. Development and Validation of a Novel Diagnostic Test for Human Brucellosis Using a Glyco-engineered Antigen Coupled to Magnetic Beads. PLoS Negl Trop Dis. 2013 Feb;7(2):e2048.

-Iwashkiw, J., Fentabil, M., Faridmoayer, A., Mills, D.C., Peppler, M. Czibener, C., Ciochini, A.E., Comerci, D.J., Ugalde, J.E. and Feldman, M.F. Exploiting the Campylobacter jejuni protein glycosylation system for glycoengineering vaccines and diagnostic tools directed against brucellosis. Microbial Cell Factories 2012 Jan 25;11:13.

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Chemtest S.A.