VETSTRIP® Brucella Bovinos

Formato dipstick

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Referencias

-Ciocchini, AE; Rey Serantes, DA; Melli, LJ; Guidolin, LS; Iwashkiw, J; Elena, S; Franco, C; Nicola, AM; Feldman, MF; Comerci, DJ and Ugalde, JE. A bacterial engineered glycoprotein as a novel antigen for diagnosis of bovine brucellosis. *Veterinary Microbiology*, Aug. 2014, 27; 172(3-4):455-65.

-Iwashkiw, J., Fentabil, M., Faridmoayer, A., Mills, D.C., Peppler, M. Czibener, C., Ciocchini, A.E., Comerci, D.J., Ugalde, J.E. and Feldman, M.F. Exploiting the *Campylobacter jejuni* protein glycosylation system for glycoengineering vaccines and diagnostic tools directed against brucellosis. *Microbial Cell Factories* 2012 Jan 25;11-13.

VETSTRIP® Brucella Bovine (serum)

Dipstick format

Kit components

<i>Strips</i>	In sealed tri-laminated pouch to protect the strips from light and moisture	20 strips (5 strips per pouch)
<i>Sample diluent</i>	Ready to use in plastic dropper bottle	1 x 7 ml
<i>Instruction manual</i>		1

Name and application

VETSTRIP® Brucella Bovine (serum)

Dipstick format

One step immunochromatographic test for the detection of specific antibodies against the O polysaccharide of the lipopolysaccharide (LPS) of Brucella in bovine serum.

Thanks to its exclusive **GlycoEng** technology, the kit has an excellent diagnostic performance, minimizes cross-reaction with other Gram-negative bacteria and allows to differentiate naturally infected animals from those vaccinated with the vaccine strain S19.

Introduction

Brucellosis is a highly contagious zoonotic disease caused by Gram-negative bacteria of the genus Brucella. The species Brucella abortus is the causative agent of brucellosis in cattle. This disease causes abortions in advanced gestation periods in pregnant females, infertility, placental retention, premature births and decreased milk production. In males, the infection can cause infertility, epididymitis and orchitis. Additionally, detection of the disease in a region or country imposes restrictions on the movement and sale of livestock. For these reasons, bovine brucellosis causes great economic losses and, therefore, should be monitored in herds to detect infection. Due to the high economic impact of the disease in animal production and health and to the risk of transmission to humans, most countries have implemented programs to control and/or eradicate brucellosis in cattle. In the absence of a vaccine against human brucellosis, disease prevention depends primarily on the control of brucellosis in animals, the natural reservoir of the disease. In particular, the control of bovine brucellosis depends primarily on vaccination with the vaccine strain S19,

in combination with the detection and slaughter of infected animals. The confirmatory diagnosis of brucellosis is performed directly by isolating the organism from blood cultures, bone marrow or other tissues. However, due to the slow growth of Brucella in primary cultures (up to 7 days), the risk involved in manipulating the bacteria and low sensitivity of its isolation, determined that the isolation of Brucella for diagnosis is not always feasible and effective. For this reason, the **laboratory diagnosis of bovine brucellosis is mainly based on the detection of specific antibodies against the infectious agent in serum and/or milk samples**. In this sense, it is essential that the serological method used allows discrimination of vaccinated from naturally infected animals, and thus avoid unnecessary slaughter of healthy animals.

Principle of the technique

The VETSTRIP® Brucella Bovine is an immunochromatographic test based on the detection of IgG antibodies directed against the O polysaccharide of the LPS of Brucella.

The test is based on the immunological capture of gold nanoparticles coated with goat anti-bovine IgG antibodies during its passage through a membrane. The IgG isotype immunoglobulins present in the serum sample react with gold particles functionalized with anti-bovine IgG antibodies and the formed complexes migrate by chromatography to the reaction zone. In this area, the antigen (purified recombinant glycoprotein) and donkey anti-goat IgG antibodies were immobilized in the test (TL) and control lines (CL), respectively. If the sample contains bovine IgG antibodies directed against the O polysaccharide of Brucella, the nanoparticle-antibody complexes are captured by the antigen immobilized on the membrane (TL) which is visualized as a pink/red band. Regardless of the presence in the sample of IgG antibod-

ies anti-Brucella O polysaccharide, the complexes that were not captured in the TL continue to migrate and are captured in the CL by the anti-goat IgG antibodies, leading also to the formation of a pink/red band. The appearance of the control line indicates that the chromatography has been developed correctly and under conditions that ensure the antigen-antibody reaction.

Kit components

-Strips in dipstick format

Reactive strips in sealed tri-laminated pouch with silica gel packets inside to protect them from light and moisture.

-Sample diluent (dropper bottle with 7 ml)

Ready to use.

-Instruction manual

Storage and expiration

Store at room temperature between 15 and 25°C.

Transport temperature: 15 to 25°C.

Keep the strips in their original packaging with the silica packets and hermetically sealed to protect them from light and moisture.

Proper storage of the kit and all reagents will allow the use until the expiration date.

Expiration: 12 months.

Do not use reagents from different lots or after the expiration date.

Materials needed but not provided

-Timer.

-Khan test tubes or flat bottom 96-well microtitre plates.

-Micropipette.

-Micropipette tips.

Precautions

1- Carefully read and strictly follow all instructions detailed in the manual included in the kit.

2- Store the kit and all reagents at room temperature (15-25°C).

3- Keep the strips in their original packaging with the silica packets and hermetically sealed to protect them from light and moisture.

4- Remove the strips from its pouch only immediately before using it.

5- Handle all materials and reagents according to the Good Laboratory Practices.

6- Serum samples should be considered potentially infectious and, therefore, the kit components that contact them should be handled with gloves and disposed as pathological waste according to legal regulations. The rest of the kit components does not represent any potential risk for health and the environment.

7- Do not use the kit beyond expiration date or if the kit components were not kept in the conditions indicated above.

8- Do not mix components or instruction manuals from different test kits batches.

Reagents and samples

-Sample diluent

Ready to use in plastic dropper bottle.

-Serum samples

Fresh, chilled or frozen serum samples, free of turbidity, can be analyzed. Samples can be stored in refrigerator for 1 or 2 days. For longer periods store at -20°C and, in these cases, the samples should be completely thawed, bring to room temperature and homogenize before performing the analysis. Samples contaminated and/or not well conserved may produce erroneous results.

Instructions for the analysis (see Figure 1)

1. Dispense 4 drops of Sample diluent (approximately 100 µl).

2. Add to the tube or plate 5 µl of serum sample (previously allow the samples equil-

ibrate to room temperature). Homogenize.

- Take one test strip from the top (violet color) and place it vertically on the tube or plate following the orientation of the arrows printed on the test strip. Close the envelope containing the remaining strips with adhesive tape.
- Wait 10 minutes, read and interpret the test result.


Figure 1: Instructions for the analysis and interpretation of results.

Interpretation of results and validity criteria of the test (see Figure 1)

Positive test result: when two pink/red bands are observed. One corresponding to the test line (TL) and the other corresponding to the control line (CL).

Note: regardless of the color intensity of the TL the result should be considered positive.

Negative test result: when only one pink/red band corresponding to the control line (CL) is observed.

Invalid result: when no bands are observed or when only one band corre-

sponding to the test line (TL) is observed. In this case, the test must be repeated in strict accordance with the procedure detailed in the kit's Instruction manual.

Any line that may appear after the 10 minutes will not have diagnostic value.

Sensitivity and specificity

To determine the diagnostic sensitivity and specificity of the assay, serum samples obtained from animals vaccinated with strain S19 and establishments free of brucellosis (negative samples), and from serologically positive animals by BPA, CELISA, FPA and CFT* (positive samples) were analyzed.

The group of negative samples included serologically negative samples for all tests (BPAT, CELISA, FPA and CFT*) as well as positive samples by BPAT but negative for the rest of the tests (CELISA, FPA and CFT).

Based on the obtained results, the diagnostic sensitivity and specificity of the assay were calculated:

Sensitivity = 97%

Specificity = 100%

*BPAT, buffered plate agglutination assay; CELISA, competitive ELISA; FPA, fluorescence polarization assay; CFT, complement fixation test.

References

-Ciocchini, AE; Rey Serantes, DA; Melli, LJ; Guidolin, LS; Iwashkiw, J; Elena, S; Franco, C; Nicola, AM; Feldman, MF; Comerci, DJ and Ugalde, JE. A bacterial engineered glycoprotein as a novel antigen for diagnosis of bovine brucellosis. *Veterinary Microbiology*, Aug. 2014, 27; 172(3-4):455-65.

-Iwashkiw, J., Fentabil, M., Faridmoayer, A., Mills, D.C., Pepler, M. Czubener, C., Ciocchini, A.E., Comerci, D.J., Ugalde, J.E. and Feldman, M.F. Exploiting the *Campylobacter jejuni* protein glycosylation system for glycoengineering vaccines and diagnostic tools directed against brucellosis. *Microbial Cell Factories* 2012 Jan 25;11:13.

00

LFIA

Bovino

 Chemtest S.A.

Av. 25 de Mayo 1021, San Martín (C.P. 1650), Buenos Aires, Argentina.
info@chemtest.net - chemtest.net