# VETSTRIP® Brucella Caninos

## Formato dipstick

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Dipstick format

Kit components

<i>Strips</i>	In sealed tri-laminated pouch to protect the strips from light and moisture	20 strips (5 strips per pouch)
<i>Sample diluent</i>	Ready to use in plastic dropper bottle	1 x 7 ml
<i>Instruction manual</i>		1

## Name and application

### VETSTRIP® Brucella Canine

Dipstick format

One step immunochromatographic test for the detection of specific antibodies against the rough lipopolysaccharide (rLPS) of *Brucella canis* in canine serum samples.

### Introduction

Brucellosis is a highly contagious zoonosis disease caused by Gram-negative bacteria of the genus *Brucella*. The species *Brucella canis* is the causative agent of brucellosis in dogs and also a human pathogen. *B. canis* causes abortions and stillbirths in pregnant bitches and epididymitis, orchitis and prostatitis in males. These clinical manifestations can lead to sterility, a factor that produces important economic losses in dog breeding kennels. Dogs that have been neutered do not have reproductive signs, but they occasionally develop other conditions such as ocular disease and discospondylitis. Among dogs the infection is transmitted through contaminated aborted fetuses, urine, milk, vaginal secretions and semen. Additionally, *B. canis* can spread through fomites. It has been demonstrated that uninfected dogs living with infected animals of the same sex can acquire the infection within 6 months. Because of its zoonotic potential, *B. canis* can also be transmitted to humans causing the disease.

The control of canine brucellosis depends primarily on the detection and separation of infected animals. The confirmatory diagnosis of brucellosis is performed directly by isolating the organism from blood and/or semen/vaginal swab culture. However, due to the slow growth of *Brucella* in primary cultures (up to 7 days), the risk involved in manipulating the bacteria and low sensitivity of its isolation, determined that the isolation of

*Brucella* for diagnosis is not always feasible and effective. For this reason, the laboratory diagnosis of canine brucellosis is mainly based on the detection of specific antibodies against the infectious agent in serum samples.

### Principle of the technique

The VETSTRIP® *Brucella Canine* is an immunochromatographic test based on the detection of IgG antibodies directed against the *Brucella canis* rough lipopolysaccharide (rLPS).

The test is based on the immunological capture of gold nanoparticles coated with goat anti-canine IgG antibodies during its passage through a membrane. The IgG isotype immunoglobulins present in the serum sample react with gold particles functionalized with anti-canine IgG antibodies and the formed complexes migrate by chromatography to the reaction zone. In this area, the antigen (LPS<sub>r</sub>) and donkey anti-goat IgG antibodies were immobilized in the test (TL) and control lines (CL), respectively.

If the sample contains canine IgG antibodies directed against *B. canis* LPS<sub>r</sub>, the nanoparticle-antibody complexes are captured by the antigen immobilized on the membrane (TL) which is visualized as a pink/red band. Regardless of the presence in the sample of IgG antibodies anti-*B. canis* LPS<sub>r</sub>, the complexes that were not captured in the TL continue to migrate and are captured in CL by the anti-goat IgG antibodies, leading also to the formation of a pink/red band. The appearance of the control line indicates that the chromatography has been developed correctly and under conditions that ensure the antigen-antibody reaction.

### Kit components

#### -Strips in dipstick format

Reactive strips in sealed tri-laminated pouch with silica gel packets inside to

protect them from light and moisture.

-Sample diluent (dropper bottle with 7 ml)

Ready to use.

-Instruction manual

### **Storage and expiration**

Store at room temperature between 15 and 25°C.

Transport temperature: 15 to 25°C.

Keep the strips in their original packaging with the silica packets and hermetically sealed to protect them from light and moisture.

Proper storage of the kit and all reagents will allow the use until the expiration date.

Expiration: 12 months.

Do not use reagents from different lots or after the expiration date.

### **Materials needed but not provided**

-Timer.

-Khan test tubes or flat bottom 96-well microtitre plates.

-Micropipette.

-Micropipette tips.

### **Precautions**

1- Carefully read and strictly follow all instructions detailed in the manual included in the kit.

2- Store the kit and all reagents at room temperature (15-25°C).

3- Keep the strips in their original packaging with the silica packets and hermetically sealed to protect them from light and moisture.

4- Remove the strips from its pouch only immediately before using it.

5- Handle all materials and reagents according to the Good Laboratory Practices.

6- Serum samples should be considered potentially infectious and, therefore, the kit components that contact them

should be handled with gloves and disposed as pathological waste according to legal regulations. The rest of the kit components does not represent any potential risk for health and the environment.

7- Do not use the kit beyond expiration date or if the kit components were not kept in the conditions indicated above.

8- Do not mix components or instruction manuals from different test kits batches.

### **Reagents and samples**

-Sample diluent

Ready to use in plastic dropper bottle.

-Serum samples

Fresh, chilled or frozen serum samples, free of turbidity, can be analyzed. Samples can be stored in refrigerator for 1 or 2 days. For longer periods store at -20°C and, in these cases, the samples should be completely thawed, bring to room temperature and homogenize before performing the analysis. Samples contaminated and/or not well conserved may produce erroneous results.

### **Instructions for the analysis (see Figure 1)**

1. Dispense 4 drops of Sample diluent (approximately 100 µl).

2. Add to the tube or plate 5 µl of serum sample (previously allow the samples equilibrate to room temperature). Homogenize.

3. Take one test strip from the top (violet color) and place it vertically on the tube or plate following the orientation of the arrows printed on the test strip. Close the envelope containing the remaining strips with adhesive tape.

4. Wait 10 minutes, read and interpret the test result.

### **Interpretation of results and validity criteria of the test (see Figure 1)**

**Positive test result:** when two pink/red bands are observed. One corresponding

to the test line (TL) and the other corresponding to the control line (CL).

Note: regardless of the color intensity of the TL the result should be considered positive.

**Negative test result:** when only one pink/red band corresponding to the control line (CL) is observed.

**Invalid result:** when no bands are observed or when only one band corresponding to the test line (TL) is observed. In this case, the test must be repeated in strict accordance with the procedure detailed in the kit's Instruction manual.

Any line that may appear after the 10 minutes will not have diagnostic value.


Figure 1: Instructions for the analysis and interpretation of results.

## Sensitivity and specificity

To validate the immunochromatographic assay VETSTRIP® Brucella Canine two different sample panels were analyzed.

**Panel I.** Samples used to determine the diagnostic sensitivity and specificity of the tests. This panel included 71 serum

samples obtained from culture-positive and serologically positive (RSAT\* and AGID\* positive) dogs ( $n = 38$  from Argentina and  $n = 33$  from Colombia), and 213 samples obtained from healthy dogs from Argentina with negative results by RSAT and AGID. The serum samples from Argentina were provided by the DILAB-SENASA and the sera from Colombia were provided by the Biogénesis group-Universidad de Antioquia.

**Panel II.** 203 encoded serum samples provided by the Programa Nacional de Tenencia Responsable y Sanidad de Perros y Gatos (ProTenencia, Ministerio de Salud de Nación, Argentina). All the samples were analyzed in parallel by RSAT (at SENASA) and the kit VETSTRIP® Brucella Canine. Samples that resulted positive by RSAT and/or the kit were additionally analyzed by AGID.

Based on the obtained results, the diagnostic sensitivity (Se) was 98.6% and the diagnostic specificity (Sp) 100%. In addition, an almost perfect correlation of VETSTRIP® Brucella Canine with RSAT (Cohen statistic,  $\kappa = 0.831$ ) was observed.

\* RSAT, rapid slide agglutination test; AGID, agar gel immunodiffusion test.

## References

-Development of Improved Enzyme-based and Lateral Flow Immunoassays for Rapid and Accurate Serodiagnosis of Canine Brucellosis. M. E. Cortina, A. Novak, L. J. Melli, S. Elena, N. Corbera, J. E. Romero, A. M. Nicola, J. E. Ugalde D. J. Comerchi and A. E. Ciochini. *Veterinary Microbiology* 2017 Sep;208:174-180.

00  
LFIA

 Chemtest S.A.  
Av. 25 de Mayo 1021, San Martín (C.P. 1650), Buenos Aires, Argentina.  
[info@chemtest.net](mailto:info@chemtest.net) - [chemtest.net](http://chemtest.net)